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Página 3

Item Name	#	ACU73																					
Auto Rerun SW																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Auto Rerun Range (Results)																							
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 30%;"></th> <th style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF Lower</th> <th style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF Higher</th> </tr> <tr> <td>Serum</td> <td style="text-align: center;">1,5</td> <td style="text-align: center;">20</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF Higher	Serum	1,5	20	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF Higher																					
Serum	1,5	20																					
Urine																							
Plasma																							
CSF																							
Dialysis																							
Other																							
Auto Rerun Condition (Absorbance)																							
Absorbance Range																							
<input type="checkbox"/> Lower <input checked="" type="radio"/> ON <input type="radio"/> OFF																							
<input type="checkbox"/> Higher <input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Prozone Range																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							

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Página 1

Item Name	#	ALB19		
DATA INFORMATION				
UNITS	g/dL			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W.Length 1	600			
Sub W.length 2	700			
METHOD	Verde de Bromocresol			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2		#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Página 3

Item Name	#	ALB19																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">1,0</td> <td style="text-align: center;">6,0</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	1,0	6,0	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
Serum	1,0	6,0																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Lower</td> <td style="width: 30%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																					
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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Página 1

Item Name	#	ALT108		
DATA INFORMATION				
UNITS	U/I			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	340			
Sub W.length 2	700			
METHOD	Cinética UV – IFCC			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>Third Mix. <input type="radio"/> OFF <input type="radio"/> ON</p> <p>R1 Blank <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1</p>	<p>FACTOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Blank correction</td> <td style="width: 30%; text-align: center;">1</td> </tr> </table>	Blank correction	1
Blank correction	1		

<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="width: 30%; text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>PROZONE CHECK</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 20%;">START</td> <td style="width: 20%;">END</td> <td style="width: 30%;">LIMIT(%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> </tr> </table> <p style="text-align: right;"> <input type="radio"/> Low <input type="radio"/> High <input type="radio"/> Low <input type="radio"/> High </p>		START	END	LIMIT(%)	FIRST				SECOND				THIRD			
0 LEVEL POINT	1																				
SPAN	3																				
	START	END	LIMIT(%)																		
FIRST																					
SECOND																					
THIRD																					

Página 3

Item Name	#	ALT108
------------------	----------	---------------

<p>Auto Rerun SW</p> <p><input type="radio"/> ON <input checked="" type="radio"/> OFF</p>	<p>Auto Rerun Condition (Absorbance)</p> <p>Absorbance Range</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Lower</td> <td style="width: 30%;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF				
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF				

<p>Auto Rerun Range (Results)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">10,0</td> <td style="text-align: center;">400</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	10,0	400	Urine			Plasma			CSF			Dialysis			Other			<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
	Lower	Higher																							
Serum	10,0	400																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

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Página 1

Item Name	#	ALT108		
DATA INFORMATION				
UNITS	U/L			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	340			
Sub W.length 2	700			
METHOD	Cinética UV – IFCC			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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ALT/GPT Liquiform

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Página 2

Item Name	#	ALT108	
ASPIRATION			
<input type="checkbox"/> KIND <input type="radio"/> Single <input checked="" type="radio"/> Double			
	VOLUME		
SAMPLE	20	µl	
REAGENT1 VOL	160		
REAGENT2	40		
DATA PROCESS			
READ			
	START	END	
MAIN	38	42	
SUB			
ABSORBANCE LIMIT			
LOW	1.100		
HIGH	2.500		
FACTOR			
Blank correction	1		
ENDPOINT LIMIT			
ENDPOINT LIMIT			
LINEAR CHECK(%)	40		
MONITOR			
0 LEVEL POINT	1		
SPAN	3		
PROZONE CHECK			
	START	END	LIMIT(%)
FIRST			
SECOND			
THIRD			
			<input type="radio"/> Low <input type="radio"/> High <input type="radio"/> Low <input type="radio"/> High

Página 3

Item Name	#	ALT108
Auto Rerun SW		
<input type="radio"/> ON <input checked="" type="radio"/> OFF		
Auto Rerun Range (Results)		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	Lower	Higher
Serum	10,0	400
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
Prozone Range		
<input type="radio"/> ON <input checked="" type="radio"/> OFF		

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Página 1

Item Name	#	AMI25		
DATA INFORMATION				
UNITS	U/I			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	405			
Sub W.length 2	700			
METHOD	CNPG3			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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AMILASE CNPG

Página 2

Item Name	#	AMI25
ASPIRATION		DATA PROCESS
KIND <input type="radio"/> Single <input type="radio"/> Double		READ
VOLUME		START END
SAMPLE	4	MAIN 22 26
REAGENT1 VOL	200	SUB
µl		
		ABSORBANCE LIMIT
		LOW -2.500
		HIGH 1.450
Third Mix. <input type="radio"/> OFF <input type="radio"/> ON		FACTOR
R1 Blank <input type="radio"/> Water Blank <input type="radio"/> R1-Blank1		Blank correction 1
		ENDPOINT LIMIT
		LINEAR CHECK(%) 40
MONITOR		PROZONE CHECK
0 LEVEL POINT 1		START END LIMIT(%)
SPAN 3		FIRST
		SECOND <input type="radio"/> Low <input type="radio"/> High
		THIRD <input type="radio"/> Low <input type="radio"/> High

Página 3

Item Name	#	AMI25
Auto Rerun SW		
<input type="radio"/> ON <input checked="" type="radio"/> OFF		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
Auto Rerun Range (Results)		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	Lower	Higher
Serum	25	2000
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF		

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Página 1

Item Name	#	AST109		
DATA INFORMATION				
UNITS	U/I			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	340			
Sub W.length 2	700			
METHOD	Cinética UV – IFCC			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>Auto Rerun Range (Results)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">10</td> <td style="text-align: center;">400</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	10	400	Urine			Plasma			CSF			Dialysis			Other			<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
	Lower	Higher																							
Serum	10	400																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

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Página 1

Item Name	#	AST109		
DATA INFORMATION				
UNITS	U/L			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	340			
Sub W.length 2	700			
METHOD	Cinética UV – IFCC			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>ASPIRATION</p> <p>KIND <input type="radio"/> Single <input checked="" type="radio"/> Double</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">VOLUME</td> <td style="width: 30%;"></td> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">20</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">160</td> </tr> <tr> <td>REAGENT2</td> <td style="text-align: center;">40</td> </tr> </table> <p style="text-align: right;">µl</p>	VOLUME		SAMPLE	20	REAGENT1 VOL	160	REAGENT2	40	<p>DATA PROCESS</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="3">READ</td> </tr> <tr> <td style="width: 30%;"></td> <td style="width: 35%;">START</td> <td style="width: 35%;">END</td> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">38</td> <td style="text-align: center;">42</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table>	READ				START	END	MAIN	38	42	SUB		
VOLUME																					
SAMPLE	20																				
REAGENT1 VOL	160																				
REAGENT2	40																				
READ																					
	START	END																			
MAIN	38	42																			
SUB																					

<p>Third Mix. <input type="radio"/> OFF <input type="radio"/> ON</p> <p>R1 Blank <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1</p>	<p>ABSORBANCE LIMIT</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">LOW</td> <td style="width: 40%; text-align: center;">1.100</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">2.500</td> </tr> </table>	LOW	1.100	HIGH	2.500
LOW	1.100				
HIGH	2.500				

<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">0 LEVEL POINT</td> <td style="width: 40%; text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>FACTOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Blank correction</td> <td style="width: 40%; text-align: center;">1</td> </tr> </table>	Blank correction	1
0 LEVEL POINT	1						
SPAN	3						
Blank correction	1						

PROZONE CHECK			
	START	END	LIMIT(%)
FIRST			
SECOND			
THIRD			

Low High
 Low High

Página 3

Item Name	#	AST109
------------------	----------	---------------

<p>Auto Rerun SW</p> <p><input type="radio"/> ON <input checked="" type="radio"/> OFF</p>	<p>Auto Rerun Condition (Absorbance)</p> <p>Absorbance Range</p> <p>Lower <input type="radio"/> ON <input checked="" type="radio"/> OFF</p> <p>Higher <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
--	---

<p>Auto Rerun Range (Results)</p> <table style="width: 100%;"> <tr> <td style="width: 25%;"></td> <td style="width: 12.5%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 12.5%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">10</td> <td style="text-align: center;">400</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	10	400	Urine			Plasma			CSF			Dialysis			Other			<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
	Lower	Higher																							
Serum	10	400																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

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Intervalo operacional: 12 mg/dL

Página 1

Item Name	#	BDLiq		
DATA INFORMATION				
UNITS	mg/dL			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W.Length 1	546			
Sub W.length 2	660			
METHOD	DCA			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="width: 30%; text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>PROZONE CHECK</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 20%; text-align: center;">START</td> <td style="width: 20%; text-align: center;">END</td> <td style="width: 30%; text-align: center;">LIMIT(%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> </tr> </table> <p style="text-align: right;"> <input type="radio"/> Low <input type="radio"/> High <input type="radio"/> Low <input type="radio"/> High </p>		START	END	LIMIT(%)	FIRST				SECOND				THIRD			
0 LEVEL POINT	1																				
SPAN	3																				
	START	END	LIMIT(%)																		
FIRST																					
SECOND																					
THIRD																					

<input type="text" value="Third Mix."/>	<input type="radio"/> OFF	<input type="radio"/> ON
<input type="text" value="R1 Blank"/>	<input checked="" type="radio"/> Water Blank	<input type="radio"/> R1-Blank1

Página 3

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="BDLiq"/>
--	--------------------------------	------------------------------------

<p>Auto Rerun SW</p> <p><input type="radio"/> ON <input checked="" type="radio"/> OFF</p>	<p>Auto Rerun Condition (Absorbance)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Absorbance Range</td> <td style="width: 30%;"></td> </tr> <tr> <td style="text-align: center;">Lower</td> <td><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td style="text-align: center;">Higher</td> <td><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Absorbance Range		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF
Absorbance Range							
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF						
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF						

<p>Auto Rerun Range (Results)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input type="radio"/> OFF</td> </tr> <tr> <td style="text-align: center;">Lower</td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">Higher</td> <td></td> <td></td> </tr> <tr> <td>Serum</td> <td></td> <td></td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF	Lower			Higher			Serum			Urine			Plasma			CSF			Dialysis			Other			<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Prozone Range</td> <td style="width: 30%;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Prozone Range	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF																												
Lower																														
Higher																														
Serum																														
Urine																														
Plasma																														
CSF																														
Dialysis																														
Other																														
Prozone Range	<input type="radio"/> ON <input checked="" type="radio"/> OFF																													

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Intervalo operacional: 30 mg/dL

Página 1

Item Name	#	BTLiq		
DATA INFORMATION				
UNITS	mg/dL			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W.Length 1	546			
Sub W.length 2				
METHOD	DCA			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Bili-T Liquiform

Página 2

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="BTLiq"/>																	
ASPIRATION																			
<input type="text" value="KIND"/> <input type="radio"/> Single <input checked="" type="radio"/> Double																			
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">SAMPLE</th> <th style="width: 50%;">VOLUME</th> </tr> <tr> <td>SAMPLE</td> <td>12.5</td> </tr> <tr> <td>REAGENT1 VOL</td> <td>200</td> </tr> <tr> <td>REAGENT2</td> <td>50</td> </tr> </table>		SAMPLE	VOLUME	SAMPLE	12.5	REAGENT1 VOL	200	REAGENT2	50	µl									
SAMPLE	VOLUME																		
SAMPLE	12.5																		
REAGENT1 VOL	200																		
REAGENT2	50																		
<input type="text" value="Third Mix."/> <input type="radio"/> OFF <input type="radio"/> ON																			
<input type="text" value="R1 Blank"/> <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1																			
MONITOR																			
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">0 LEVEL POINT</td> <td style="width: 50%;">1</td> </tr> <tr> <td>SPAN</td> <td>3</td> </tr> </table>		0 LEVEL POINT	1	SPAN	3														
0 LEVEL POINT	1																		
SPAN	3																		
DATA PROCESS																			
READ																			
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 30%;">START</th> <th style="width: 70%;">END</th> </tr> <tr> <td>MAIN</td> <td>53 54</td> </tr> <tr> <td>SUB</td> <td>30 31</td> </tr> </table>		START	END	MAIN	53 54	SUB	30 31												
START	END																		
MAIN	53 54																		
SUB	30 31																		
ABSORBANCE LIMIT																			
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">LOW</td> <td style="width: 50%;">0</td> </tr> <tr> <td>HIGH</td> <td>3.000</td> </tr> </table>		LOW	0	HIGH	3.000														
LOW	0																		
HIGH	3.000																		
FACTOR																			
<input type="text" value="Blank correction"/> <input type="text" value="0.8095"/>																			
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">ENDPOINT LIMIT</td> <td style="width: 50%;">0.010</td> </tr> <tr> <td>LINEAR CHECK(%)</td> <td>0</td> </tr> </table>		ENDPOINT LIMIT	0.010	LINEAR CHECK(%)	0														
ENDPOINT LIMIT	0.010																		
LINEAR CHECK(%)	0																		
PROZONE CHECK																			
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 25%;"></th> <th style="width: 25%;">START</th> <th style="width: 25%;">END</th> <th style="width: 25%;">LIMIT(%)</th> </tr> </thead> <tbody> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> </tr> </tbody> </table>					START	END	LIMIT(%)	FIRST				SECOND				THIRD			
	START	END	LIMIT(%)																
FIRST																			
SECOND																			
THIRD																			
		<input type="radio"/> Low	<input type="radio"/> High																
		<input type="radio"/> Low	<input type="radio"/> High																

Página 3

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="BTLiq"/>																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%;"> <tr> <td></td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;"><input type="text" value="Lower"/></td> <td style="text-align: center;"><input type="text" value="Higher"/></td> </tr> <tr> <td><input type="text" value="Serum"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="Urine"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="Plasma"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="CSF"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="Dialysis"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="Other"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		<input type="text" value="Lower"/>	<input type="text" value="Higher"/>	<input type="text" value="Serum"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="Urine"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="Plasma"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="CSF"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="Dialysis"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="Other"/>	<input type="text"/>	<input type="text"/>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	<input type="text" value="Lower"/>	<input type="text" value="Higher"/>																								
<input type="text" value="Serum"/>	<input type="text"/>	<input type="text"/>																								
<input type="text" value="Urine"/>	<input type="text"/>	<input type="text"/>																								
<input type="text" value="Plasma"/>	<input type="text"/>	<input type="text"/>																								
<input type="text" value="CSF"/>	<input type="text"/>	<input type="text"/>																								
<input type="text" value="Dialysis"/>	<input type="text"/>	<input type="text"/>																								
<input type="text" value="Other"/>	<input type="text"/>	<input type="text"/>																								
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<input type="text" value="Lower"/>		<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
<input type="text" value="Higher"/>		<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
Prozone Range																										
<input type="text" value="Prozone Range"/>		<input type="radio"/> ON <input checked="" type="radio"/> OFF																								

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Página 1

Item Name	#	Ca95		
DATA INFORMATION				
UNITS	mg/dL			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W.Length 1	660			
Sub W.length 2				
METHOD	Arsenazo			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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LINEAR CHECK(%)

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				<input type="radio"/> Low <input type="radio"/> High
SECOND				<input type="radio"/> Low <input type="radio"/> High
THIRD				<input type="radio"/> Low <input type="radio"/> High

Página 3

Item Name	#	Ca95
------------------	----------	-------------

Auto Rerun SW

ON OFF

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower ON OFF

Higher ON OFF

Auto Rerun Range (Results)

	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	Lower	Higher
Serum	3	17
Urine		
Plasma		
CSF		
Dialysis		
Other		

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Página 1

Item Name	#	Ca90		
DATA INFORMATION				
UNITS	mg/dl			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W.Length 1	570			
Sub W.length 2	660			
METHOD	CPC			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Página 3

Item Name	#	Ca90																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 33%;"></td> <td style="width: 33%; text-align: center;"> <input type="radio"/> ON <input checked="" type="radio"/> OFF </td> <td style="width: 33%; text-align: center;"> <input type="radio"/> ON <input checked="" type="radio"/> OFF </td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td><input type="checkbox"/> Serum</td> <td style="text-align: center;">3</td> <td style="text-align: center;">16</td> </tr> <tr> <td><input type="checkbox"/> Urine</td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Plasma</td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> CSF</td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Dialysis</td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	<input type="checkbox"/> Serum	3	16	<input type="checkbox"/> Urine			<input type="checkbox"/> Plasma			<input type="checkbox"/> CSF			<input type="checkbox"/> Dialysis			<input type="checkbox"/> Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
<input type="checkbox"/> Serum	3	16																								
<input type="checkbox"/> Urine																										
<input type="checkbox"/> Plasma																										
<input type="checkbox"/> CSF																										
<input type="checkbox"/> Dialysis																										
<input type="checkbox"/> Other																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<input type="checkbox"/> Lower <input type="radio"/> ON <input checked="" type="radio"/> OFF																										
<input type="checkbox"/> Higher <input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Prozone Range																										
<input type="checkbox"/> ON <input checked="" type="radio"/> OFF																										

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<p>ASPIRATION</p> <p>KIND <input type="radio"/> Single <input type="radio"/> Double</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">VOLUME</td> <td style="width: 30%;"></td> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">10</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">200</td> </tr> </table> <p style="text-align: right;">µl</p>	VOLUME		SAMPLE	10	REAGENT1 VOL	200	<p>DATA PROCESS</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="3">READ</td> </tr> <tr> <td style="width: 30%;"></td> <td style="width: 35%;">START</td> <td style="width: 35%;">END</td> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">35</td> <td style="text-align: center;">54</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">ABSORBANCE LIMIT</td> </tr> <tr> <td style="width: 60%;">LOW</td> <td style="text-align: center;">-0.010</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">0.375</td> </tr> </table>	READ				START	END	MAIN	35	54	SUB			ABSORBANCE LIMIT		LOW	-0.010	HIGH	0.375
VOLUME																									
SAMPLE	10																								
REAGENT1 VOL	200																								
READ																									
	START	END																							
MAIN	35	54																							
SUB																									
ABSORBANCE LIMIT																									
LOW	-0.010																								
HIGH	0.375																								

<p>Third Mix. <input type="radio"/> OFF <input type="radio"/> ON</p> <p>R1 Blank <input type="radio"/> Water Blank <input type="radio"/> R1-Blank1</p>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">FACTOR</td> </tr> <tr> <td style="width: 60%;">Blank correction</td> <td style="text-align: center;">1</td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">ENDPOINT LIMIT</td> </tr> <tr> <td style="width: 60%;">LINEAR CHECK(%)</td> <td style="text-align: center;">40</td> </tr> </table>	FACTOR		Blank correction	1	ENDPOINT LIMIT		LINEAR CHECK(%)	40
FACTOR									
Blank correction	1								
ENDPOINT LIMIT									
LINEAR CHECK(%)	40								

<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>PROZONE CHECK</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 25%;"></td> <td style="width: 25%;">START</td> <td style="width: 25%;">END</td> <td style="width: 25%;">LIMIT(%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> </tr> </table> <p style="text-align: right;"> <input type="radio"/> Low <input type="radio"/> High <input type="radio"/> Low <input type="radio"/> High </p>		START	END	LIMIT(%)	FIRST				SECOND				THIRD			
0 LEVEL POINT	1																				
SPAN	3																				
	START	END	LIMIT(%)																		
FIRST																					
SECOND																					
THIRD																					

Página 3

Item Name	#	MB78
-----------	---	------

<p>Auto Rerun SW</p> <p><input type="radio"/> ON <input checked="" type="radio"/> OFF</p>	<p>Auto Rerun Condition (Absorbance)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">Absorbance Range</td> </tr> <tr> <td style="width: 60%;">Lower</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Absorbance Range		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF
Absorbance Range							
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF						
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF						

<p>Auto Rerun Range (Results)</p> <table style="width: 100%;"> <tr> <td style="width: 25%;"></td> <td style="width: 25%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 25%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 25%;"></td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> <td></td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">0</td> <td style="text-align: center;">500</td> <td></td> </tr> <tr> <td>Urine</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF			Lower	Higher		Serum	0	500		Urine				Plasma				CSF				Dialysis				Other				<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Prozone Range</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Prozone Range	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																																	
	Lower	Higher																																	
Serum	0	500																																	
Urine																																			
Plasma																																			
CSF																																			
Dialysis																																			
Other																																			
Prozone Range	<input type="radio"/> ON <input checked="" type="radio"/> OFF																																		

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Página 1

Item Name	#	CK77		
DATA INFORMATION				
UNITS	U/I			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	340			
Sub W.length 2	700			
METHOD				
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2		#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>ASPIRATION</p> <p>KIND <input type="radio"/> Single <input type="radio"/> Double</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">VOLUME</td> <td style="width: 30%;"></td> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">4</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">200</td> </tr> </table> <p style="text-align: right;">µl</p>	VOLUME		SAMPLE	4	REAGENT1 VOL	200	<p>DATA PROCESS</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="3">READ</td> </tr> <tr> <td style="width: 30%;"></td> <td style="width: 35%;">START</td> <td style="width: 35%;">END</td> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">20</td> <td style="text-align: center;">28</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">ABSORBANCE LIMIT</td> </tr> <tr> <td style="width: 60%;">LOW</td> <td style="text-align: center;">-0.010</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">0.900</td> </tr> </table>	READ				START	END	MAIN	20	28	SUB			ABSORBANCE LIMIT		LOW	-0.010	HIGH	0.900
VOLUME																									
SAMPLE	4																								
REAGENT1 VOL	200																								
READ																									
	START	END																							
MAIN	20	28																							
SUB																									
ABSORBANCE LIMIT																									
LOW	-0.010																								
HIGH	0.900																								

<p>Third Mix. <input type="radio"/> OFF <input type="radio"/> ON</p> <p>R1 Blank <input type="radio"/> Water Blank <input type="radio"/> R1-Blank1</p>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">FACTOR</td> </tr> <tr> <td style="width: 60%;">Blank correction</td> <td style="text-align: center;">1</td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">ENDPOINT LIMIT</td> </tr> <tr> <td style="width: 60%;">LINEAR CHECK(%)</td> <td style="text-align: center;">40</td> </tr> </table>	FACTOR		Blank correction	1	ENDPOINT LIMIT		LINEAR CHECK(%)	40
FACTOR									
Blank correction	1								
ENDPOINT LIMIT									
LINEAR CHECK(%)	40								

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<p>Auto Rerun Range (Results)</p> <table style="width: 100%;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">25</td> <td style="text-align: center;">2000</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	25	2000	Urine			Plasma			CSF			Dialysis			Other			<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
	Lower	Higher																							
Serum	25	2000																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

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Página 3

Item Name	#	HDL13																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">10</td> <td style="text-align: center;">200</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	10	200	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
Serum	10	200																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Lower</td> <td style="width: 30%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																					
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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Página 1

Item Name	#	CRE96		
DATA INFORMATION				
UNITS	mg/dl			
DECIMALS	2			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	505			
Sub W.length 2	570			
METHOD	Labtest			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER - 0.25		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>PROZONE CHECK</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td></td> <td style="width: 15%;">START</td> <td style="width: 15%;">END</td> <td style="width: 15%;">LIMIT(%)</td> <td></td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> <td><input type="radio"/> Low <input type="radio"/> High</td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> <td><input type="radio"/> Low <input type="radio"/> High</td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> <td><input type="radio"/> Low <input type="radio"/> High</td> </tr> </table>		START	END	LIMIT(%)		FIRST				<input type="radio"/> Low <input type="radio"/> High	SECOND				<input type="radio"/> Low <input type="radio"/> High	THIRD				<input type="radio"/> Low <input type="radio"/> High
0 LEVEL POINT	1																								
SPAN	3																								
	START	END	LIMIT(%)																						
FIRST				<input type="radio"/> Low <input type="radio"/> High																					
SECOND				<input type="radio"/> Low <input type="radio"/> High																					
THIRD				<input type="radio"/> Low <input type="radio"/> High																					

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<p>Auto Rerun Range (Results)</p> <table style="width: 100%;"> <tr> <td></td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">0,2</td> <td style="text-align: center;">12</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	0,2	12	Urine			Plasma			CSF			Dialysis			Other			<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
	Lower	Higher																							
Serum	0,2	12																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

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<p>ASPIRATION</p> <p>KIND <input type="radio"/> Single <input checked="" type="radio"/> Double</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">VOLUME</td> <td style="width: 30%;"></td> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">10</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">160</td> </tr> <tr> <td>REAGENT2 VOL</td> <td style="text-align: center;">40</td> </tr> </table> <p style="text-align: right;">µl</p>	VOLUME		SAMPLE	10	REAGENT1 VOL	160	REAGENT2 VOL	40	<p>DATA PROCESS</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="3">READ</td> </tr> <tr> <td style="width: 30%;"></td> <td style="width: 35%;">START</td> <td style="width: 35%;">END</td> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">36</td> <td style="text-align: center;">48</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table>	READ				START	END	MAIN	36	48	SUB		
VOLUME																					
SAMPLE	10																				
REAGENT1 VOL	160																				
REAGENT2 VOL	40																				
READ																					
	START	END																			
MAIN	36	48																			
SUB																					

<p>Third Mix. <input type="radio"/> OFF <input type="radio"/> ON</p> <p>R1 Blank <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1</p>	<p>ABSORBANCE LIMIT</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">LOW</td> <td style="text-align: center;">-2.500</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">0.600</td> </tr> </table>	LOW	-2.500	HIGH	0.600
LOW	-2.500				
HIGH	0.600				

<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">0 LEVEL POINT</td> <td style="text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>FACTOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Blank correction</td> <td style="text-align: center;">1</td> </tr> </table>	Blank correction	1
0 LEVEL POINT	1						
SPAN	3						
Blank correction	1						

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<p>Auto Rerun SW</p> <p><input type="radio"/> ON <input checked="" type="radio"/> OFF</p>	<p>Auto Rerun Condition (Absorbance)</p> <p>Absorbance Range</p> <p style="text-align: center;">Lower <input type="radio"/> ON <input checked="" type="radio"/> OFF</p> <p style="text-align: center;">Higher <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
--	---

<p>Auto Rerun Range (Results)</p> <table style="width: 100%;"> <tr> <td style="width: 25%;"></td> <td style="width: 12.5%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 12.5%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 50%;"></td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> <td></td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">0,2</td> <td style="text-align: center;">12</td> <td></td> </tr> <tr> <td>Urine</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF			Lower	Higher		Serum	0,2	12		Urine				Plasma				CSF				Dialysis				Other				<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																															
	Lower	Higher																															
Serum	0,2	12																															
Urine																																	
Plasma																																	
CSF																																	
Dialysis																																	
Other																																	

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Página 1

Item Name	#	FeLiq		
DATA INFORMATION				
UNITS	µg/dL			
DECIMALS	0			
ANALYSIS				
TYPE	END			
Main W.Length 1	570			
Sub W.length 2				
METHOD	Ferrozine			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2		#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Fe Liquiform

Página 2

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="FeLiq"/>																										
ASPIRATION <input type="text" value="KIND"/> <input type="radio"/> Single <input checked="" type="radio"/> Double <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;"></td> <td style="width: 30%; text-align: center;">VOLUME</td> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">25</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">200</td> </tr> <tr> <td>REAGENT2</td> <td style="text-align: center;">50</td> </tr> </table> <div style="text-align: right; margin-top: 5px;">µl</div>			VOLUME	SAMPLE	25	REAGENT1 VOL	200	REAGENT2	50	DATA PROCESS <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="3">READ</td> </tr> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;">START</td> <td style="width: 35%; text-align: center;">END</td> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">53</td> <td style="text-align: center;">54</td> </tr> <tr> <td>SUB</td> <td style="text-align: center;">30</td> <td style="text-align: center;">31</td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <tr> <td colspan="2">ABSORBANCE LIMIT</td> </tr> <tr> <td style="width: 60%;">LOW</td> <td style="text-align: center;">-0.010</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">0.550</td> </tr> </table>	READ				START	END	MAIN	53	54	SUB	30	31	ABSORBANCE LIMIT		LOW	-0.010	HIGH	0.550
	VOLUME																											
SAMPLE	25																											
REAGENT1 VOL	200																											
REAGENT2	50																											
READ																												
	START	END																										
MAIN	53	54																										
SUB	30	31																										
ABSORBANCE LIMIT																												
LOW	-0.010																											
HIGH	0.550																											
<input type="text" value="Third Mix."/> <input type="radio"/> OFF <input type="radio"/> ON <input type="text" value="R1 Blank"/> <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">FACTOR</td> </tr> <tr> <td style="width: 60%;">Blank correction</td> <td style="text-align: center;">0.8182</td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <tr> <td style="width: 60%;">ENDPOINT LIMIT</td> <td style="text-align: center;">0.010</td> </tr> <tr> <td>LINEAR CHECK(%)</td> <td></td> </tr> </table>	FACTOR		Blank correction	0.8182	ENDPOINT LIMIT	0.010	LINEAR CHECK(%)																			
FACTOR																												
Blank correction	0.8182																											
ENDPOINT LIMIT	0.010																											
LINEAR CHECK(%)																												
MONITOR <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">0 LEVEL POINT</td> <td style="text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>		0 LEVEL POINT	1	SPAN	3	PROZONE CHECK <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 20%; text-align: center;">START</td> <td style="width: 20%; text-align: center;">END</td> <td style="width: 30%; text-align: center;">LIMIT(%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> </tr> </table> <div style="margin-top: 5px;"> <input type="radio"/> Low <input type="radio"/> High <input checked="" type="radio"/> Low <input type="radio"/> High </div>		START	END	LIMIT(%)	FIRST				SECOND				THIRD									
0 LEVEL POINT	1																											
SPAN	3																											
	START	END	LIMIT(%)																									
FIRST																												
SECOND																												
THIRD																												

Página 3

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="FeLiq"/>																										
Auto Rerun SW <input type="radio"/> ON <input checked="" type="radio"/> OFF		Auto Rerun Condition (Absorbance) <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">Absorbance Range</td> </tr> <tr> <td style="width: 60%; text-align: center;">Lower</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td style="text-align: center;">Higher</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Absorbance Range		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																				
Absorbance Range																												
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																											
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																											
Auto Rerun Range (Results) <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input type="radio"/> OFF</td> </tr> <tr> <td style="text-align: center;">Lower</td> <td></td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td></td> <td style="text-align: center;">1000</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>			<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF	Lower		Higher	Serum		1000	Urine			Plasma			CSF			Dialysis			Other			<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Prozone Range</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Prozone Range	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF																										
Lower		Higher																										
Serum		1000																										
Urine																												
Plasma																												
CSF																												
Dialysis																												
Other																												
Prozone Range	<input type="radio"/> ON <input checked="" type="radio"/> OFF																											

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Página 1

Item Name	#	FAL79		
DATA INFORMATION				
UNITS	U/I			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	405			
Sub W.length 2	700			
METHOD				
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Página 3

Item Name	#	FAL79																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<input type="checkbox"/> Lower		<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
<input type="checkbox"/> Higher		<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
Auto Rerun Range (Results)																										
<table style="width: 100%; border-collapse: collapse;"> <tr> <td></td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">20</td> <td style="text-align: center;">1500</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	20	1500	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
Serum	20	1500																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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Página 1

Item Name	#	FOS12		
DATA INFORMATION				
UNITS	mg/dl			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W.Length 1	340			
Sub W.length 2	660			
METHOD				
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>ASPIRATION</p> <p>KIND <input type="radio"/> Single <input type="radio"/> Double</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">VOLUME</td> <td style="width: 30%;"></td> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">3</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">300</td> </tr> </table> <p style="text-align: right;">µl</p>	VOLUME		SAMPLE	3	REAGENT1 VOL	300	<p>DATA PROCESS</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="3">READ</td> </tr> <tr> <td style="width: 30%;"></td> <td style="width: 35%;">START</td> <td style="width: 35%;">END</td> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">25</td> <td style="text-align: center;">26</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table>	READ				START	END	MAIN	25	26	SUB		
VOLUME																			
SAMPLE	3																		
REAGENT1 VOL	300																		
READ																			
	START	END																	
MAIN	25	26																	
SUB																			

<p>Third Mix. <input type="radio"/> OFF <input type="radio"/> ON</p> <p>R1 Blank <input type="radio"/> Water Blank <input type="radio"/> R1-Blank1</p>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">ABSORBANCE LIMIT</td> </tr> <tr> <td style="width: 60%;">LOW</td> <td style="text-align: center;">0.600</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">2.300</td> </tr> </table>	ABSORBANCE LIMIT		LOW	0.600	HIGH	2.300
ABSORBANCE LIMIT							
LOW	0.600						
HIGH	2.300						

<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">0 LEVEL POINT</td> <td style="text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>FACTOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Blank correction</td> <td style="text-align: center;">1</td> </tr> </table>	Blank correction	1
0 LEVEL POINT	1						
SPAN	3						
Blank correction	1						

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<p>Auto Rerun SW</p> <p><input type="radio"/> ON <input checked="" type="radio"/> OFF</p>	<p>Auto Rerun Condition (Absorbance)</p> <p>Absorbance Range</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Lower</td> <td style="width: 40%;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF				
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF				

<p>Auto Rerun Range (Results)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 25%;"></td> <td style="width: 12.5%;"><input type="radio"/> ON <input type="radio"/> OFF</td> <td style="width: 12.5%;"></td> <td style="width: 12.5%;"><input type="radio"/> ON <input type="radio"/> OFF</td> </tr> <tr> <td></td> <td>Lower</td> <td>Higher</td> <td></td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">2,5</td> <td style="text-align: center;">20</td> <td></td> </tr> <tr> <td>Urine</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input type="radio"/> OFF		<input type="radio"/> ON <input type="radio"/> OFF		Lower	Higher		Serum	2,5	20		Urine				Plasma				CSF				Dialysis				Other				<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input type="radio"/> OFF		<input type="radio"/> ON <input type="radio"/> OFF																														
	Lower	Higher																															
Serum	2,5	20																															
Urine																																	
Plasma																																	
CSF																																	
Dialysis																																	
Other																																	

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Página 3

Item Name	#	FRU97																					
Auto Rerun SW																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Auto Rerun Range (Results)																							
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"> <input type="radio"/> ON <input type="radio"/> OFF Lower </td> <td style="width: 35%; text-align: center;"> <input type="radio"/> ON <input type="radio"/> OFF Higher </td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">20</td> <td style="text-align: center;">800</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input type="radio"/> OFF Lower	<input type="radio"/> ON <input type="radio"/> OFF Higher	Serum	20	800	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input type="radio"/> OFF Lower	<input type="radio"/> ON <input type="radio"/> OFF Higher																					
Serum	20	800																					
Urine																							
Plasma																							
CSF																							
Dialysis																							
Other																							
Auto Rerun Condition (Absorbance)																							
Absorbance Range																							
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%; text-align: center;">Lower</td> <td style="width: 30%;"> <input type="radio"/> ON <input checked="" type="radio"/> OFF </td> </tr> <tr> <td style="text-align: center;">Higher</td> <td> <input type="radio"/> ON <input checked="" type="radio"/> OFF </td> </tr> </table>			Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																	
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																						
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																						
Prozone Range																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							

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<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>PROZONE CHECK</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 20%;">START</td> <td style="width: 20%;">END</td> <td style="width: 30%;">LIMIT(%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> </tr> </table> <p style="text-align: right;"> <input type="radio"/> Low <input type="radio"/> High <input type="radio"/> Low <input type="radio"/> High </p>		START	END	LIMIT(%)	FIRST				SECOND				THIRD			
0 LEVEL POINT	1																				
SPAN	3																				
	START	END	LIMIT(%)																		
FIRST																					
SECOND																					
THIRD																					

Item Name	#	GGT105
------------------	----------	---------------

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<p>Auto Rerun Range (Results)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">5</td> <td style="text-align: center;">700</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	5	700	Urine			Plasma			CSF			Dialysis			Other			<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
	Lower	Higher																							
Serum	5	700																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

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Página 1

Item Name	#	GHK85		
DATA INFORMATION				
UNITS	mg/dl			
DECIMALS	0			
ANALYSIS				
TYPE	END			
Main W.Length 1	340			
Sub W.length 2	380			
METHOD				
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE		Linear		
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Página 3

Item Name	#	GHK85																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%; border-collapse: collapse;"> <tr> <td></td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">50</td> <td style="text-align: center;">700</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	50	700	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
Serum	50	700																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<input type="checkbox"/> Lower <input type="radio"/> ON <input checked="" type="radio"/> OFF																										
<input type="checkbox"/> Higher <input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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Página 1

Item Name	#	GLI84		
DATA INFORMATION				
UNITS	mg/dl			
DECIMALS	0			
ANALYSIS				
TYPE	END			
Main W.Length 1	505			
Sub W.length 2	660			
METHOD	GOD – Trinder			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Página 3

Item Name	#	GLI84																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">50</td> <td style="text-align: center;">500</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF		Lower	Higher	Serum	50	500	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF																								
	Lower	Higher																								
Serum	50	500																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Lower</td> <td style="width: 30%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																					
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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Página 3

Item Name	#	HDL98																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%;"> <tr> <td></td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">30</td> <td style="text-align: center;">200</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	30	200	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
Serum	30	200																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<input type="checkbox"/> Lower <input type="radio"/> ON <input checked="" type="radio"/> OFF																										
<input type="checkbox"/> Higher <input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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Página 1

Item Name	#	IBC		
DATA INFORMATION				
UNITS	µg/dL			
DECIMALS	0			
ANALYSIS				
TYPE	END			
Main W.Length 1	570			
Sub W.length 2				
METHOD	Labtest Ferrozine			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2		#5		
#3		#6		
NORMAL RANGE				
	MALE	FEMALE		
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p><input type="text" value="Third Mix."/> <input type="radio"/> OFF <input type="radio"/> ON</p> <p><input type="text" value="R1 Blank"/> <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1</p>	<p>ABSORBANCE LIMIT</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">LOW</td> <td style="width: 40%; text-align: center;">0.350</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">0.750</td> </tr> </table>	LOW	0.350	HIGH	0.750
LOW	0.350				
HIGH	0.750				

<p>MONITOR</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="width: 30%; text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>	0 LEVEL POINT	1	SPAN	3	<p>FACTOR</p> <p><input type="text" value="Blank correction"/> <input type="text" value="0.8065"/></p>
0 LEVEL POINT	1				
SPAN	3				

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<p>Auto Rerun SW</p> <p><input type="radio"/> ON <input checked="" type="radio"/> OFF</p>	<p>Auto Rerun Condition (Absorbance)</p> <p>Absorbance Range</p> <p><input type="text" value="Lower"/> <input type="radio"/> ON <input checked="" type="radio"/> OFF</p> <p><input type="text" value="Higher"/> <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
--	---

<p>Auto Rerun Range (Results)</p> <table style="width: 100%;"> <tr> <td></td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td style="text-align: center;">Lower</td> <td><input type="text"/></td> <td style="text-align: center;">Higher</td> </tr> <tr> <td><input type="text" value="Serum"/></td> <td><input type="text"/></td> <td style="text-align: center;">500</td> </tr> <tr> <td><input type="text" value="Urine"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="Plasma"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="CSF"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="Dialysis"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td><input type="text" value="Other"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Lower	<input type="text"/>	Higher	<input type="text" value="Serum"/>	<input type="text"/>	500	<input type="text" value="Urine"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="Plasma"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="CSF"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="Dialysis"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="Other"/>	<input type="text"/>	<input type="text"/>	<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Lower	<input type="text"/>	Higher																							
<input type="text" value="Serum"/>	<input type="text"/>	500																							
<input type="text" value="Urine"/>	<input type="text"/>	<input type="text"/>																							
<input type="text" value="Plasma"/>	<input type="text"/>	<input type="text"/>																							
<input type="text" value="CSF"/>	<input type="text"/>	<input type="text"/>																							
<input type="text" value="Dialysis"/>	<input type="text"/>	<input type="text"/>																							
<input type="text" value="Other"/>	<input type="text"/>	<input type="text"/>																							

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Página 1

Item Name	#	LDH86		
DATA INFORMATION				
UNITS	U/I			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	340			
Sub W.length 2	700			
METHOD	Piruvato-Lactato			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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LDH Liquiform

Página 2

Item Name	#	LDH86																				
ASPIRATION																						
<input type="checkbox"/> KIND <input type="radio"/> Single <input checked="" type="radio"/> Double																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 70%;">VOLUME</th> <th style="width: 30%;"></th> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">4</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">160</td> </tr> <tr> <td>REAGENT2</td> <td style="text-align: center;">40</td> </tr> </table>		VOLUME		SAMPLE	4	REAGENT1 VOL	160	REAGENT2	40	µl												
VOLUME																						
SAMPLE	4																					
REAGENT1 VOL	160																					
REAGENT2	40																					
<input type="checkbox"/> Third Mix. <input type="radio"/> OFF <input type="radio"/> ON																						
<input type="checkbox"/> R1 Blank <input type="radio"/> Water Blank <input type="radio"/> R1-Blank1																						
MONITOR																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>		0 LEVEL POINT	1	SPAN	3																	
0 LEVEL POINT	1																					
SPAN	3																					
DATA PROCESS																						
READ																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 30%;">START</th> <th style="width: 70%;">END</th> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">37 41</td> </tr> <tr> <td>SUB</td> <td></td> </tr> </table>		START	END	MAIN	37 41	SUB																
START	END																					
MAIN	37 41																					
SUB																						
ABSORBANCE LIMIT																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">LOW</td> <td style="text-align: center;">1.1</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">2.5</td> </tr> </table>		LOW	1.1	HIGH	2.5																	
LOW	1.1																					
HIGH	2.5																					
FACTOR																						
<input type="checkbox"/> Blank correction <input type="checkbox"/> 1																						
ENDPOINT LIMIT																						
<input type="checkbox"/> ENDPOINT LIMIT		40																				
PROZONE CHECK																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;"></th> <th style="width: 15%;">START</th> <th style="width: 15%;">END</th> <th style="width: 30%;">LIMIT(%)</th> <th style="width: 10%;"></th> </tr> </thead> <tbody> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> <td rowspan="2" style="text-align: center;"><input type="radio"/> Low <input type="radio"/> High</td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> <td style="text-align: center;"><input type="radio"/> Low <input type="radio"/> High</td> </tr> </tbody> </table>					START	END	LIMIT(%)		FIRST				<input type="radio"/> Low <input type="radio"/> High	SECOND				THIRD				<input type="radio"/> Low <input type="radio"/> High
	START	END	LIMIT(%)																			
FIRST				<input type="radio"/> Low <input type="radio"/> High																		
SECOND																						
THIRD				<input type="radio"/> Low <input type="radio"/> High																		

Página 3

Item Name	#	LDH86																					
Auto Rerun SW																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Auto Rerun Condition (Absorbance)																							
Absorbance Range																							
<input type="checkbox"/> Lower <input type="radio"/> ON <input checked="" type="radio"/> OFF																							
<input type="checkbox"/> Higher <input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Auto Rerun Range (Results)																							
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;"></th> <th style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF Lower</th> <th style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF Higher</th> </tr> </thead> <tbody> <tr> <td>Serum</td> <td style="text-align: center;">200</td> <td style="text-align: center;">2000</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </tbody> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF Higher	Serum	200	2000	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF Higher																					
Serum	200	2000																					
Urine																							
Plasma																							
CSF																							
Dialysis																							
Other																							
Prozone Range																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							

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<input radio"="" type="text" value="Third Mix.</td> <td><input type="/> OFF <input type="radio"/> ON	
<input type="text" value="R1 Blank"/>	<input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1

MONITOR	
0 LEVEL POINT	1
SPAN	3

FACTOR			
<input type="text" value="Blank correction"/>	0.7526		
ENDPOINT LIMIT		0.020	
LINEAR CHECK(%)			

PROZONE CHECK			
	START	END	LIMIT(%)
FIRST			
SECOND			
THIRD			
			<input type="radio"/> Low <input type="radio"/> High
			<input type="radio"/> Low <input type="radio"/> High

Página 3

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="LDL"/>
--	--------------------------------	----------------------------------

Auto Rerun SW	Auto Rerun Condition (Absorbance)						
<input type="radio"/> ON <input checked="" type="radio"/> OFF	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2" style="text-align: center;">Absorbance Range</td> </tr> <tr> <td style="width: 60%;"><input type="text" value="Lower"/></td> <td style="text-align: right;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td><input type="text" value="Higher"/></td> <td style="text-align: right;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>	Absorbance Range		<input type="text" value="Lower"/>	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="text" value="Higher"/>	<input type="radio"/> ON <input checked="" type="radio"/> OFF
Absorbance Range							
<input type="text" value="Lower"/>	<input type="radio"/> ON <input checked="" type="radio"/> OFF						
<input type="text" value="Higher"/>	<input type="radio"/> ON <input checked="" type="radio"/> OFF						

Auto Rerun Range (Results)	Prozone Range																								
<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td></td> <td style="width: 20%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 20%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td style="text-align: center;">Lower</td> <td></td> <td style="text-align: center;">Higher</td> </tr> <tr> <td style="text-align: center;">Serum</td> <td style="text-align: center;">6.6</td> <td style="text-align: center;">992</td> </tr> <tr> <td style="text-align: center;">Urine</td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">Plasma</td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">CSF</td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">Dialysis</td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Lower		Higher	Serum	6.6	992	Urine			Plasma			CSF			Dialysis			Other			<input type="radio"/> ON <input checked="" type="radio"/> OFF
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Lower		Higher																							
Serum	6.6	992																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

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NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Lipase Liquiform

Página 2

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="LIP"/>	
ASPIRATION			
<input type="text" value="KIND"/> <input type="radio"/> Single <input checked="" type="radio"/> Double			
<input type="text" value="VOLUME"/>			
<input type="text" value="SAMPLE"/>	<input type="text" value="2.5"/>	µl	
<input type="text" value="REAGENT1 VOL"/>	<input type="text" value="155"/>	µl	
<input type="text" value="REAGENT2"/>	<input type="text" value="95"/>	µl	
DATA PROCESS			
READ			
	<input type="text" value="START"/>	<input type="text" value="END"/>	
<input type="text" value="MAIN"/>	<input type="text" value="50"/>	<input type="text" value="51"/>	
<input type="text" value="SUB"/>	<input type="text" value="38"/>	<input type="text" value="39"/>	
ABSORBANCE LIMIT			
<input type="text" value="LOW"/>	<input type="text" value="-0.010"/>		
<input type="text" value="HIGH"/>	<input type="text" value="2.500"/>		
FACTOR			
<input type="text" value="Blank correction"/>	<input type="text" value="1"/>		
<input type="text" value="ENDPOINT LIMIT"/>		<input type="text" value="0.020"/>	
<input type="text" value="LINEAR CHECK(%)"/>			
PROZONE CHECK			
	<input type="text" value="START"/>	<input type="text" value="END"/>	<input type="text" value="LIMIT(%)"/>
<input type="text" value="FIRST"/>			
<input type="text" value="SECOND"/>			
<input type="text" value="THIRD"/>			
			<input type="radio"/> Low <input type="radio"/> High
			<input type="radio"/> Low <input type="radio"/> High
MONITOR			
<input type="text" value="0 LEVEL POINT"/>	<input type="text" value="1"/>		
<input type="text" value="SPAN"/>	<input type="text" value="3"/>		
Third Mix. <input type="radio"/> OFF <input type="radio"/> ON			
R1 Blank <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1			

Página 3

<input type="text" value="Item Name"/>	<input type="text" value="#"/>	<input type="text" value="LDL"/>
Auto Rerun SW		
<input type="radio"/> ON <input checked="" type="radio"/> OFF		
Auto Rerun Range (Results)		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	<input type="text" value="Lower"/>	<input type="text" value="Higher"/>
<input type="text" value="Serum"/>	<input type="text" value="3.0"/>	<input type="text" value="300"/>
<input type="text" value="Urine"/>	<input type="text"/>	<input type="text"/>
<input type="text" value="Plasma"/>	<input type="text"/>	<input type="text"/>
<input type="text" value="CSF"/>	<input type="text"/>	<input type="text"/>
<input type="text" value="Dialysis"/>	<input type="text"/>	<input type="text"/>
<input type="text" value="Other"/>	<input type="text"/>	<input type="text"/>
Auto Rerun Condition (Absorbance)		
Absorbance Range		
<input type="text" value="Lower"/>	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
<input type="text" value="Higher"/>	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF		

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Página 1

Item Name	#	Mg50		
DATA INFORMATION				
UNITS	mg/dl			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W.Length 1	505			
Sub W.length 2				
METHOD	Labtest			
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Página 3

Item Name	#	Mg50																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Range (Results)																										
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">2</td> <td style="text-align: center;">4.5</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	2	4.5	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
Serum	2	4.5																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Lower</td> <td style="width: 30%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																					
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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PROTEÍNAS TOTAIS

Página 2

Item Name	#	PRO99	
ASPIRATION			
KIND	<input type="radio"/> Single <input type="radio"/> Double		
SAMPLE	VOLUME		
REAGENT1 VOL	200	µl	
DATA PROCESS			
READ			
	START	END	
MAIN	53	54	
SUB			
ABSORBANCE LIMIT			
LOW	-0.300		
HIGH	0.800		
FACTOR			
Blank correction	1		
ENDPOINT LIMIT			
ENDPOINT LIMIT	0.010		
LINEAR CHECK(%)			
MONITOR			
0 LEVEL POINT	1		
SPAN	3		
PROZONE CHECK			
	START	END	LIMIT(%)
FIRST			
SECOND			
THIRD			
			<input type="radio"/> Low <input type="radio"/> High
			<input type="radio"/> Low <input type="radio"/> High

Página 3

Item Name	#	PRO99
Auto Rerun SW		
<input type="radio"/> ON <input checked="" type="radio"/> OFF		
Auto Rerun Range (Results)		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF
	Lower	Higher
Serum	1	14
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF	
Prozone Range		
<input type="radio"/> ON <input checked="" type="radio"/> OFF		

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SENSIPROT

Página 2

Item Name	#	SEN36																				
ASPIRATION																						
<input type="checkbox"/> KIND <input checked="" type="radio"/> Single <input type="radio"/> Double																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">VOLUME</td> <td style="width: 30%;"></td> </tr> <tr> <td>SAMPLE</td> <td style="text-align: center;">10</td> </tr> <tr> <td>REAGENT1 VOL</td> <td style="text-align: center;">200</td> </tr> </table>		VOLUME		SAMPLE	10	REAGENT1 VOL	200	µl														
VOLUME																						
SAMPLE	10																					
REAGENT1 VOL	200																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Third Mix.</td> <td style="width: 30%;"> <input checked="" type="radio"/> OFF <input type="radio"/> ON </td> </tr> <tr> <td>R1 Blank</td> <td> <input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1 </td> </tr> </table>			Third Mix.	<input checked="" type="radio"/> OFF <input type="radio"/> ON	R1 Blank	<input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1																
Third Mix.	<input checked="" type="radio"/> OFF <input type="radio"/> ON																					
R1 Blank	<input checked="" type="radio"/> Water Blank <input type="radio"/> R1-Blank1																					
MONITOR																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">0 LEVEL POINT</td> <td style="width: 30%; text-align: center;">1</td> </tr> <tr> <td>SPAN</td> <td style="text-align: center;">3</td> </tr> </table>		0 LEVEL POINT	1	SPAN	3																	
0 LEVEL POINT	1																					
SPAN	3																					
DATA PROCESS																						
READ																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;">START</td> <td style="width: 70%;">END</td> </tr> <tr> <td>MAIN</td> <td style="text-align: center;">53 54</td> </tr> <tr> <td>SUB</td> <td></td> </tr> </table>		START	END	MAIN	53 54	SUB																
START	END																					
MAIN	53 54																					
SUB																						
ABSORBANCE LIMIT																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">LOW</td> <td style="width: 30%; text-align: center;">0.090</td> </tr> <tr> <td>HIGH</td> <td style="text-align: center;">1.000</td> </tr> </table>		LOW	0.090	HIGH	1.000																	
LOW	0.090																					
HIGH	1.000																					
FACTOR																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Blank correction</td> <td style="width: 30%; text-align: center;">1</td> </tr> </table>		Blank correction	1																			
Blank correction	1																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">ENDPOINT LIMIT</td> <td style="width: 30%; text-align: center;">0.010</td> </tr> <tr> <td>LINEAR CHECK(%)</td> <td></td> </tr> </table>		ENDPOINT LIMIT	0.010	LINEAR CHECK(%)																		
ENDPOINT LIMIT	0.010																					
LINEAR CHECK(%)																						
PROZONE CHECK																						
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;"></th> <th style="width: 15%;">START</th> <th style="width: 15%;">END</th> <th style="width: 30%;">LIMIT(%)</th> <th style="width: 10%;"></th> </tr> </thead> <tbody> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> <td rowspan="2" style="text-align: center;"> <input type="radio"/> Low <input type="radio"/> High </td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td></td> <td style="text-align: center;"> <input checked="" type="radio"/> Low <input type="radio"/> High </td> </tr> </tbody> </table>					START	END	LIMIT(%)		FIRST				<input type="radio"/> Low <input type="radio"/> High	SECOND				THIRD				<input checked="" type="radio"/> Low <input type="radio"/> High
	START	END	LIMIT(%)																			
FIRST				<input type="radio"/> Low <input type="radio"/> High																		
SECOND																						
THIRD				<input checked="" type="radio"/> Low <input type="radio"/> High																		

Página 3

Item Name	#	SEN36																					
Auto Rerun SW																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
Auto Rerun Range (Results)																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;"></th> <th style="width: 35%; text-align: center;"> <input type="radio"/> ON <input checked="" type="radio"/> OFF Lower </th> <th style="width: 35%; text-align: center;"> <input type="radio"/> ON <input checked="" type="radio"/> OFF Higher </th> </tr> </thead> <tbody> <tr> <td>Serum</td> <td></td> <td></td> </tr> <tr> <td>Urine</td> <td style="text-align: center;">1</td> <td style="text-align: center;">100</td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td style="text-align: center;">15</td> <td style="text-align: center;">100</td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </tbody> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF Higher	Serum			Urine	1	100	Plasma			CSF	15	100	Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF Higher																					
Serum																							
Urine	1	100																					
Plasma																							
CSF	15	100																					
Dialysis																							
Other																							
Auto Rerun Condition (Absorbance)																							
Absorbance Range																							
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Lower</td> <td style="width: 30%; text-align: center;"> <input type="radio"/> ON <input checked="" type="radio"/> OFF </td> </tr> <tr> <td>Higher</td> <td style="text-align: center;"> <input type="radio"/> ON <input checked="" type="radio"/> OFF </td> </tr> </table>		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																		
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																						
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																						
Prozone Range																							
<input type="radio"/> ON <input checked="" type="radio"/> OFF																							

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Página 3

Item Name	#	TGR87																								
Auto Rerun SW																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										
Auto Rerun Condition (Absorbance)																										
Absorbance Range																										
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 70%;">Lower</td> <td style="width: 30%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td>Higher</td> <td style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> </table>		Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																					
Lower	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Higher	<input type="radio"/> ON <input checked="" type="radio"/> OFF																									
Auto Rerun Range (Results)																										
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</th> <th style="text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</th> </tr> <tr> <th></th> <th style="text-align: center;">Lower</th> <th style="text-align: center;">Higher</th> </tr> </thead> <tbody> <tr> <td>Serum</td> <td style="text-align: center;">10</td> <td style="text-align: center;">1100</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </tbody> </table>				<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	10	1100	Urine			Plasma			CSF			Dialysis			Other		
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																								
	Lower	Higher																								
Serum	10	1100																								
Urine																										
Plasma																										
CSF																										
Dialysis																										
Other																										
Prozone Range																										
<input type="radio"/> ON <input checked="" type="radio"/> OFF																										

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Página 1

Item Name	#	UR104		
DATA INFORMATION				
UNITS	mg/dl			
DECIMALS	0			
ANALYSIS				
TYPE	RATE			
Main W.Length 1	340			
Sub W.length 2	700			
METHOD				
CORRELAÇÃO				
Y =	SLOPE 1	X + INTER 0		
CALIBRATION				
TYPE	Linear			
STANDARD				
#1	@	#4		
#2	@	#5		
#3		#6		
NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

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<p>Auto Rerun Range (Results)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;"></td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> <td style="width: 35%; text-align: center;"><input type="radio"/> ON <input checked="" type="radio"/> OFF</td> </tr> <tr> <td></td> <td style="text-align: center;">Lower</td> <td style="text-align: center;">Higher</td> </tr> <tr> <td>Serum</td> <td style="text-align: center;">15</td> <td style="text-align: center;">300</td> </tr> <tr> <td>Urine</td> <td></td> <td></td> </tr> <tr> <td>Plasma</td> <td></td> <td></td> </tr> <tr> <td>CSF</td> <td></td> <td></td> </tr> <tr> <td>Dialysis</td> <td></td> <td></td> </tr> <tr> <td>Other</td> <td></td> <td></td> </tr> </table>		<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF		Lower	Higher	Serum	15	300	Urine			Plasma			CSF			Dialysis			Other			<p>Prozone Range <input type="radio"/> ON <input checked="" type="radio"/> OFF</p>
	<input type="radio"/> ON <input checked="" type="radio"/> OFF	<input type="radio"/> ON <input checked="" type="radio"/> OFF																							
	Lower	Higher																							
Serum	15	300																							
Urine																									
Plasma																									
CSF																									
Dialysis																									
Other																									

Página 3