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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="ÁCIDO ÚRICO"/>	Test Code: <input type="text" value="ACU73"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/> ▼				
Reporting Unit:	<input type="text" value="mg/dL"/>	Decimal Points	<input type="text" value="2"/> ▲ ▼		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended				
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		16 ▲ ▼ ■ 33 ▲ ▼			
		Photometric <input type="text" value="2 Wavelength"/> ▼			
		Wavelength		Primary	<input type="text" value="510"/> ▼
				Secondary	<input type="text" value="660"/> ▼

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ÁCIDO ÚRICO Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="2.5"/></td> <td><input type="text" value="7"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="1.5"/></td> <td><input type="text" value="6"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="2.5"/></td> <td><input type="text" value="7"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="2.5"/>	<input type="text" value="7"/>	Female:	<input type="text" value="1.5"/>	<input type="text" value="6"/>	Other:	<input type="text" value="2.5"/>	<input type="text" value="7"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="2.5"/>	<input type="text" value="7"/>														
Female:	<input type="text" value="1.5"/>	<input type="text" value="6"/>														
Other:	<input type="text" value="2.5"/>	<input type="text" value="7"/>														
Valid Range: <input type="text" value="0"/> <input type="text" value="20"/>		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>												
Serum Index Limit Hemolysis: <input type="text" value="0"/> Icterus: <input type="text" value="0"/> Lipemia: <input type="text" value="0"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: <input type="text" value="M-Point"/> ▼ Repts: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days													
Calibrator Selection <table border="1"> <thead> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td><input type="text" value="2"/> ▲▼</td> <td><input type="text" value="1"/> # <input type="text" value="Calibra 1H"/></td> <td><input type="text" value="@1"/></td> </tr> <tr> <td><input type="text" value="2"/></td> <td><input type="text" value="2"/> # <input type="text" value="Calibra 2H"/></td> <td><input type="text" value="@2"/></td> </tr> </tbody> </table>			N°	Calibrator	Conc	<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # <input type="text" value="Calibra 1H"/>	<input type="text" value="@1"/>	<input type="text" value="2"/>	<input type="text" value="2"/> # <input type="text" value="Calibra 2H"/>	<input type="text" value="@2"/>	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Reps Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc											
<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # <input type="text" value="Calibra 1H"/>	<input type="text" value="@1"/>											
<input type="text" value="2"/>	<input type="text" value="2"/> # <input type="text" value="Calibra 2H"/>	<input type="text" value="@2"/>											
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle										

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="ALBUMINA"/>	Test Code: <input type="text" value="ALB19"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology		
Reporting Unit:	<input type="text" value="g/dL"/>	Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:		
Decimal Points	<input type="text" value="2"/>	16 <input type="text"/> <input type="text" value="0"/>		
Photometric	<input type="text" value="2 Wavelength"/>			
Wavelength	Primary	<input type="text" value="600"/>		
	Secondary	<input type="text" value="700"/>		

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ALBUMINA

Measure	Sample	Reagent	Ranges	Calibration																																																																											
<p>Sample Volumes</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">Sample Vol.</th> <th style="text-align: center;">Dilution Sample Vol.</th> <th style="text-align: center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">2 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> <td style="text-align: center;">0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/></td> </tr> </tbody> </table> <p>Diluent Code <input style="width: 100px;" type="text" value="Water"/></p> <p>Diluent Warming Limit <input style="width: 50px;" type="text" value="0"/> <input type="button" value="▲"/> <input type="button" value="▼"/> Tests</p>			Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	2 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	2	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	3	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	4	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	<p>Condition</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">1</th> <th style="text-align: center;">2</th> <th style="text-align: center;">3</th> <th style="text-align: center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> </tbody> </table>				1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	2 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>																																																																												
2	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>																																																																												
3	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>																																																																												
4	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 ALB	200 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	15 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/> tests	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/> days
R2	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/> tests	0 <input type="text"/> <input type="button" value="▲"/> <input type="button" value="▼"/> days
		<input checked="" type="checkbox"/> Stirring		
		<input type="checkbox"/> Stirring		

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ÁLBUMINA

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td>3.5</td> <td>5.5</td> </tr> <tr> <td>Female:</td> <td>3.5</td> <td>5.5</td> </tr> <tr> <td>Other:</td> <td>3.5</td> <td>5.5</td> </tr> </table>			Lower	Upper	Male:	3.5	5.5	Female:	3.5	5.5	Other:	3.5	5.5	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	3.5	5.5														
Female:	3.5	5.5														
Other:	3.5	5.5														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td>0</td> </tr> <tr> <td>Icterus</td> <td>0</td> </tr> <tr> <td>Lipemia</td> <td>0</td> </tr> </table>		Hemolysis	0	Icterus	0	Lipemia	0	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 3500 mAbs		Correction Constant Slope: 1 Intercept: 0						
Hemolysis	0															
Icterus	0															
Lipemia	0															
Valid Range: 0 6		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 0 Equation: none ▼ Judge Point: 0 ▲▼		Non-Linear Limit: 0 %												
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off																

Measure	Sample	Reagent	Ranges	Calibration											
Calibration Method: M-Point ▼ Repeats: 3 ▲▼ Curve Type: Linear ▼ Stability: 0 ▲▼ days															
Calibrator Selection <table border="1"> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> <tr> <td>2 ▲▼</td> <td>1 #</td> <td>Calibra 2H</td> <td>@₁</td> </tr> <tr> <td>2 #</td> <td>2 ▲▼</td> <td>Calibra 1H</td> <td>@₂</td> </tr> </table>			N°	Calibrator	Conc	2 ▲▼	1 #	Calibra 2H	@ ₁	2 #	2 ▲▼	Calibra 1H	@ ₂	R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Reps Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit: 0 %	
N°	Calibrator	Conc													
2 ▲▼	1 #	Calibra 2H	@ ₁												
2 #	2 ▲▼	Calibra 1H	@ ₂												
<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle															

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="ALT/GPT"/>	Test Code: <input type="text" value="ALT108"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology			
Reporting Unit:	<input type="text" value="U/L"/>	Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate			
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:			
	Decimal Points: <input type="text" value="0"/>	21 <input type="text"/> 25 <input type="text"/>			
		Photometric: <input type="text" value="2 Wavelength"/>			
		Wavelength: Primary <input type="text" value="340"/> Secondary <input type="text" value="700"/>			

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ALT/GPT Liquiform

Measure	Sample	Reagent	Ranges	Calibration																																																																											
<p>Sample Volumes</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">Sample Vol.</th> <th style="text-align: center;">Dilution Sample Vol.</th> <th style="text-align: center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">18</td> <td style="text-align: center;">0</td> <td style="text-align: center;">0</td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">0</td> <td style="text-align: center;">0</td> <td style="text-align: center;">0</td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">0</td> <td style="text-align: center;">0</td> <td style="text-align: center;">0</td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">0</td> <td style="text-align: center;">0</td> <td style="text-align: center;">0</td> </tr> </tbody> </table> <p style="margin-top: 10px;">Diluent Code Water</p> <p>Diluent Warming Limit 0 Tests</p>			Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	18	0	0	2	0	0	0	3	0	0	0	4	0	0	0	<p>Condition</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">1</th> <th style="text-align: center;">2</th> <th style="text-align: center;">3</th> <th style="text-align: center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> </tbody> </table>				1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	18	0	0																																																																												
2	0	0	0																																																																												
3	0	0	0																																																																												
4	0	0	0																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 ALT1	145	0	15 tests	0 days
R2 ALT2	35	0	15 tests	0 days

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ALT/GPT Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="11"/></td> <td><input type="text" value="45"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="10"/></td> <td><input type="text" value="37"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="11"/></td> <td><input type="text" value="45"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="11"/>	<input type="text" value="45"/>	Female:	<input type="text" value="10"/>	<input type="text" value="37"/>	Other:	<input type="text" value="11"/>	<input type="text" value="45"/>	Reaction Slope <input checked="" type="radio"/> Negative <input type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="11"/>	<input type="text" value="45"/>														
Female:	<input type="text" value="10"/>	<input type="text" value="37"/>														
Other:	<input type="text" value="11"/>	<input type="text" value="45"/>														
Valid Range: <input type="text" value="3.5"/> <input type="text" value="400"/>		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>												
Serum Index Limit Hemolysis: <input type="text" value="0"/> Icterus: <input type="text" value="0"/> Lipemia: <input type="text" value="0"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="40"/> % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: <input type="text" value="M-Point"/> ▼ Repts: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days													
Calibrator Selection <table border="1"> <thead> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td><input type="text" value="2"/> ▲▼</td> <td><input type="text" value="1"/> # <input type="text" value="Calibra 1H"/></td> <td><input type="text" value="@1"/></td> </tr> <tr> <td><input type="text" value="2"/></td> <td><input type="text" value="2"/> # <input type="text" value="Calibra 2H"/></td> <td><input type="text" value="@2"/></td> </tr> </tbody> </table>			N°	Calibrator	Conc	<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # <input type="text" value="Calibra 1H"/>	<input type="text" value="@1"/>	<input type="text" value="2"/>	<input type="text" value="2"/> # <input type="text" value="Calibra 2H"/>	<input type="text" value="@2"/>	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="-100"/> mAbs Cal Reps Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc											
<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # <input type="text" value="Calibra 1H"/>	<input type="text" value="@1"/>											
<input type="text" value="2"/>	<input type="text" value="2"/> # <input type="text" value="Calibra 2H"/>	<input type="text" value="@2"/>											
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle										

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="AMILASE"/>	Test Code: <input type="text" value="AMI25"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="U/L"/>	Decimal Points			
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="text" value="0"/>			
	<input type="radio"/> Extended	<input type="button" value="▲"/>			
		<input type="button" value="▼"/>			
		Methodology			
		Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="9"/>	<input type="text" value="13"/>		
		<input type="button" value="▲"/>	<input type="button" value="▲"/>		
		<input type="button" value="▼"/>	<input type="button" value="▼"/>		
		Photometric	<input type="text" value="2 Wavelength"/>		
		Wavelength	Primary	<input type="text" value="405"/>	<input type="button" value="▼"/>
			Secondary	<input type="text" value="700"/>	<input type="button" value="▼"/>

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 AMI	200 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/> tests	0 <input type="text"/> days
		<input checked="" type="checkbox"/> Stirring		
		<input type="checkbox"/> Stirring		

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="AST/GOT"/>	Test Code: <input type="text" value="AST109"/>		
		<input type="button" value="Print"/> <input type="button" value="Cancel"/>		
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>			
Reporting Unit:	<input type="text" value="U/L"/>	Decimal Points		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	<input type="text" value="0"/>		
		Methodology		
		Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate		
		Measuring Point:		
		<input type="text" value="21"/>	<input type="text" value="25"/>	
		Photometric <input type="text" value="2 Wavelength"/>		
		Wavelength	Primary	<input type="text" value="340"/>
			Secondary	<input type="text" value="700"/>

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 AST1 <input type="text"/>	145 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2 AST2 <input type="text"/>	35 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
		■ Stirring		

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="BILI-D"/>	Test Code: <input type="text" value="BD93"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology		
Reporting Unit:	<input type="text" value="mg/dL"/>	Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:		
Decimal Points	<input type="text" value="2"/>	16 <input type="text"/> <input type="text" value="33"/>		
		Photometric <input type="text" value="2 Wavelength"/>		
		Wavelength Primary <input type="text" value="546"/>		
		Secondary <input type="text" value="660"/>		

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 BILI-D1	180 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2 BILI-D2	45 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
		■ Stirring		

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Bili-D Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td>0</td> <td>0.4</td> </tr> <tr> <td>Female:</td> <td>0</td> <td>0.4</td> </tr> <tr> <td>Other:</td> <td>0</td> <td>0.4</td> </tr> </table>			Lower	Upper	Male:	0	0.4	Female:	0	0.4	Other:	0	0.4	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	0	0.4														
Female:	0	0.4														
Other:	0	0.4														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td>0</td> </tr> <tr> <td>Icterus</td> <td>0</td> </tr> <tr> <td>Lipemia</td> <td>0</td> </tr> </table>		Hemolysis	0	Icterus	0	Lipemia	0	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 3500 mAbs		Correction Constant Slope: 1 Intercept: 0						
Hemolysis	0															
Icterus	0															
Lipemia	0															
Valid Range: 0 12		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 0 Equation: none Judge Point: 0		Non-Linear Limit: 0 % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: M-Point Repts: 3 Curve Type: Linear Stability: 0 days													
Calibrator Selection <table border="1"> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> <tr> <td>2</td> <td>Calibra 1H</td> <td>@1</td> </tr> <tr> <td>2</td> <td>Calibra 2H</td> <td>@2</td> </tr> </table>			N°	Calibrator	Conc	2	Calibra 1H	@1	2	Calibra 2H	@2	R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Reps Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit: 0 %	
N°	Calibrator	Conc											
2	Calibra 1H	@1											
2	Calibra 2H	@2											
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle										

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="Ca Arsenazo"/>	Test Code: <input type="text" value="Ca95"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="mg/dL"/>	Decimal Points	<input type="text" value="2"/>		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended				
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		24 <input type="button" value="▲"/> <input type="button" value="▼"/> <input type="checkbox"/> 0 <input type="button" value="▲"/> <input type="button" value="▼"/>			
		Photometric <input type="text" value="1 Wavelength"/>			
		Wavelength Primary <input type="text" value="660"/>			
		Secondary <input type="text"/>			

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Measure	Sample	Reagent	Ranges	Calibration	
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability	
R1 Ca	200 <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="button" value="▲"/> <input type="button" value="▼"/>	<input checked="" type="checkbox"/> Stirring	15 <input type="button" value="▲"/> <input type="button" value="▼"/> tests	0 <input type="button" value="▲"/> <input type="button" value="▼"/> days
R2	<input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="button" value="▲"/> <input type="button" value="▼"/>	<input type="checkbox"/> Stirring	0 <input type="button" value="▲"/> <input type="button" value="▼"/> tests	0 <input type="button" value="▲"/> <input type="button" value="▼"/> days

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Ca Arsenazo Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td>8.8</td> <td>11</td> </tr> <tr> <td>Female:</td> <td>8.8</td> <td>11</td> </tr> <tr> <td>Other:</td> <td>8.8</td> <td>11</td> </tr> </table>			Lower	Upper	Male:	8.8	11	Female:	8.8	11	Other:	8.8	11	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	8.8	11														
Female:	8.8	11														
Other:	8.8	11														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td>0</td> </tr> <tr> <td>Icterus</td> <td>0</td> </tr> <tr> <td>Lipemia</td> <td>0</td> </tr> </table>		Hemolysis	0	Icterus	0	Lipemia	0	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 3500 mAbs		Correction Constant Slope: 1 Intercept: 0						
Hemolysis	0															
Icterus	0															
Lipemia	0															
Valid Range: 0 17		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 0 Equation: none Judge Point: 0		Non-Linear Limit: 0 %												
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off																

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: M-Point Repeats: 3 Curve Type: Linear Stability: 0 days													
Calibrator Selection <table border="1"> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> <tr> <td>2</td> <td>Calibra 2H</td> <td>@1</td> </tr> <tr> <td>2</td> <td>Calibra 1H</td> <td>@2</td> </tr> </table>			N°	Calibrator	Conc	2	Calibra 2H	@1	2	Calibra 1H	@2	R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Repts Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit: 0 %	
N°	Calibrator	Conc											
2	Calibra 2H	@1											
2	Calibra 1H	@2											
<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle													

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CÁLCIO Liquiform

Measure	Sample	Reagent	Ranges	Calibration																																																																											
Sample Volumes <table style="width:100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align:center;">Sample Vol.</th> <th style="text-align:center;">Dilution Sample Vol.</th> <th style="text-align:center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align:center;">1</td> <td style="text-align:center;">4 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align:center;">2</td> <td style="text-align:center;">0 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align:center;">3</td> <td style="text-align:center;">0 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align:center;">4</td> <td style="text-align:center;">0 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> <td style="text-align:center;">0 <input type="text"/></td> </tr> </tbody> </table> Diluent Code <input type="text" value="Water"/> Diluent Warming Limit <input type="text" value="0"/> Tests			Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	4 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	3	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	4	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	Condition <table style="width:100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align:center;">1</th> <th style="text-align:center;">2</th> <th style="text-align:center;">3</th> <th style="text-align:center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align:center;">●</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align:center;">●</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> <td style="text-align:center;">○</td> </tr> </tbody> </table>				1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	4 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
3	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
4	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 CaL1	150 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2 CaL2	50 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
		■ Stirring		

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CÁLCIO Liquiform

Measure	Sample	Reagent	Ranges	Calibration
Normal Range Lower Upper Male: 8.8 11 Female: 8.8 11 Other: 8.8 11 Valid Range: 0 16		Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
Serum Index Limit Hemolysis 0 Icterus 0 Lipemia 0		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit 3500 mAbs		Correction Constant Slope: 1 Intercept 0
		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit 0 Equation none ▼ Judge Point 0 ▲▼		Non-Linear Limit: 0 % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off

Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: M-Point ▼ Repls: 3 ▲▼ Curve Type: Linear ▼ Stability 0 ▲▼ days				
Calibrator Selection N° Calibrator Conc 2 1 # Calibra 2H @ ₁ 2 # Calibra 1H @ ₂			R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Repls Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit 0 %	
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle	

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="CK-MB"/>	Test Code: <input type="text" value="CKMB78"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/> ▼	Methodology			
Reporting Unit:	<input type="text" value="U/L"/>	Decimal Points	<input type="radio"/> End Point <input checked="" type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	<input type="text" value="1"/> ▲ ▼	Measuring Point:		
			<input type="text" value="25"/> ▲ ▼	<input type="text" value="33"/> ▲ ▼	
			Photometric	<input type="text" value="2 Wavelength"/> ▼	
			Wavelength	Primary	<input type="text" value="340"/> ▼
				Secondary	<input type="text" value="700"/> ▼

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CK MB Liquiform

Measure	Sample	Reagent	Ranges	Calibration																																																																											
<p>Sample Volumes</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">Sample Vol.</th> <th style="text-align: center;">Dilution Sample Vol.</th> <th style="text-align: center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">9 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> </tbody> </table> <p>Diluent Code <input type="text" value="Water"/></p> <p>Diluent Warming Limit <input type="text" value="0"/> <input type="text" value="▲"/> Tests <input type="text" value="▼"/></p>			Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	9 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	2	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	3	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	4	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	<p>Condition</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">1</th> <th style="text-align: center;">2</th> <th style="text-align: center;">3</th> <th style="text-align: center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> </tbody> </table>				1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	9 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
2	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
3	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
4	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 CKMB1 <input type="text" value="▲"/> <input type="text" value="▼"/>	145 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	<input checked="" type="checkbox"/> Stirring 15 <input type="text" value="▲"/> <input type="text" value="▼"/> tests	0 <input type="text" value="▲"/> <input type="text" value="▼"/> days
R2 CKMB2 <input type="text" value="▲"/> <input type="text" value="▼"/>	35 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	<input checked="" type="checkbox"/> Stirring 15 <input type="text" value="▲"/> <input type="text" value="▼"/> tests	0 <input type="text" value="▲"/> <input type="text" value="▼"/> days

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CK MB Liquiform

Measure	Sample	Reagent	Ranges	Calibration
Normal Range Lower Upper Male: 0 24 Female: 0 24 Other: 0 24 Valid Range: 0 500		Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
Serum Index Limit Hemolysis 0 Icterus 0 Lipemia 0		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit 3500 mAbs		Correction Constant Slope: 1 Intercept 0
		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit 0 Equation none ▼ Judge Point 0 ▲▼		Non-Linear Limit: 40 % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off

Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: 1- Point ▼ Repeats: 3 ▲▼ Curve Type: Linear ▼ Stability: 0 ▲▼ days				
Calibrator Selection N° 1 1 # 1 Calibrator Conc CAL CK @		R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Repts Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit 0 %		
		<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle		

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="CK TOTAL"/>	Test Code: <input type="text" value="CK77"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>			
Reporting Unit:	<input type="text" value="U/L"/>	Decimal Points	<input type="text" value="0"/>	
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="radio"/> Extended		
Methodology				
Type:		<input type="radio"/> End Point <input checked="" type="radio"/> Rate		
Measuring Point:				
25		32		
Photometric		<input type="text" value="2 Wavelength"/>		
Wavelength		Primary	<input type="text" value="340"/>	
		Secondary	<input type="text" value="700"/>	

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CK NAC Liquiform

Measure	Sample	Reagent	Ranges	Calibration
Sample Volumes				
	Sample Vol.	Dilution Sample Vol.	Diluent Vol.	
1	4 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	
2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	
3	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	
4	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	
Diluent Code		Water		
Diluent Warming Limit		0 <input type="text"/>	Tests	
Condition				
			1	2
			3	4
First Run			●	○
Below N-Range			○	○
Above N-Range			○	○
Panic L			○	○
Panic H			○	○
Noise			○	○
Prozone			○	○
HIGH!			○	○
ABS!			○	○
Samp. Vol Reduction			●	○

Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 CK1	160 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2 CK2	40 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
		■ Stirring		

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="COLESTEROL"/>	Test Code: <input type="text" value="COL76"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology		
Reporting Unit:	<input type="text" value="mg/dL"/>	Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:		
Decimal Points	<input type="text" value="0"/>	20 <input type="text"/>		
	<input type="button" value="▲"/> <input type="button" value="▼"/>	Photometric <input type="text" value="2 Wavelength"/>		
		Wavelength		
		Primary	<input type="text" value="510"/>	<input type="button" value="▼"/>
		Secondary	<input type="text" value="660"/>	<input type="button" value="▼"/>

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COLESTEROL Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table style="width:100%; border-collapse: collapse;"> <tr> <td style="width:50%;"></td> <td style="text-align:center;">Lower</td> <td style="text-align:center;">Upper</td> </tr> <tr> <td>Male:</td> <td style="text-align:center;"><input type="text" value="10"/></td> <td style="text-align:center;"><input type="text" value="200"/></td> </tr> <tr> <td>Female:</td> <td style="text-align:center;"><input type="text" value="10"/></td> <td style="text-align:center;"><input type="text" value="200"/></td> </tr> <tr> <td>Other:</td> <td style="text-align:center;"><input type="text" value="10"/></td> <td style="text-align:center;"><input type="text" value="200"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="10"/>	<input type="text" value="200"/>	Female:	<input type="text" value="10"/>	<input type="text" value="200"/>	Other:	<input type="text" value="10"/>	<input type="text" value="200"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="10"/>	<input type="text" value="200"/>														
Female:	<input type="text" value="10"/>	<input type="text" value="200"/>														
Other:	<input type="text" value="10"/>	<input type="text" value="200"/>														
Serum Index Limit <table style="width:100%; border-collapse: collapse;"> <tr> <td style="width:50%;">Hemolysis</td> <td style="text-align:center;"><input type="text" value="0"/></td> </tr> <tr> <td>Icterus</td> <td style="text-align:center;"><input type="text" value="0"/></td> </tr> <tr> <td>Lipemia</td> <td style="text-align:center;"><input type="text" value="0"/></td> </tr> </table>		Hemolysis	<input type="text" value="0"/>	Icterus	<input type="text" value="0"/>	Lipemia	<input type="text" value="0"/>	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>						
Hemolysis	<input type="text" value="0"/>															
Icterus	<input type="text" value="0"/>															
Lipemia	<input type="text" value="0"/>															
Valid Range: <input type="text" value="0"/> <input type="text" value="500"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> %												
				Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration																				
Calibration Method: <input type="text" value="M-Point"/> ▼ Repeats: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days																								
Calibrator Selection <table style="width:100%; border-collapse: collapse;"> <tr> <td style="width:10%;">N°</td> <td style="width:10%;"></td> <td style="width:40%;">Calibrator</td> <td style="width:10%;">Conc</td> </tr> <tr> <td style="text-align:center;">2</td> <td style="text-align:center;">1</td> <td style="text-align:center;">#</td> <td style="text-align:center;">#</td> </tr> <tr> <td></td> <td></td> <td style="text-align:center;">Calibra 2H</td> <td style="text-align:center;">@₁</td> </tr> <tr> <td></td> <td style="text-align:center;">2</td> <td style="text-align:center;">#</td> <td style="text-align:center;">#</td> </tr> <tr> <td></td> <td></td> <td style="text-align:center;">Calibra 1H</td> <td style="text-align:center;">@₂</td> </tr> </table>			N°		Calibrator	Conc	2	1	#	#			Calibra 2H	@ ₁		2	#	#			Calibra 1H	@ ₂	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Reps Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°		Calibrator	Conc																					
2	1	#	#																					
		Calibra 2H	@ ₁																					
	2	#	#																					
		Calibra 1H	@ ₂																					
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle																					

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="CREATININA"/>	Test Code: <input type="text" value="CRE96"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology		
Reporting Unit:	<input type="text" value="mg/dL"/>	Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:		
Decimal Points	<input type="text" value="2"/>	<input type="text" value="4"/>	<input type="text" value="11"/>	
				Photometric <input type="text" value="2 Wavelength"/>
				Wavelength Primary <input type="text" value="510"/>
				Secondary <input type="text" value="570"/>

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 CRE	180	0	15 tests	0 days
			<input checked="" type="checkbox"/> Stirring	
R2	0	0	0 tests	0 days
			<input type="checkbox"/> Stirring	

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CREATININA K

Measure	Sample	Reagent	Ranges	Calibration
Normal Range Lower Upper Male: 0.7 1.2 Female: 0.53 1 Other: 0.7 1.2 Valid Range: 0.2 12		Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
Serum Index Limit Hemolysis 0 Icterus 0 Lipemia 0		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit 3500 mAbs		Correction Constant Slope: 1 Intercept - 0.25
		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit 0 Equation none ▼ Judge Point 0 ▲▼		Non-Linear Limit: 40 % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off

Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: M-Point ▼ Repls: 3 ▲▼ Curve Type: Linear ▼ Stability 0 ▲▼ days				
Calibrator Selection N° Calibrator Conc 2 1 # Calibra 1H @1 2 2 # Calibra 2H @2			R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Repls Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit 0 %	
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle	

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="Fe Liquiform"/>	Test Code: <input type="text" value="Fe91"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/> ▼				
Reporting Unit:	<input type="text" value="µg/dL"/>	Decimal Points	<input type="text" value="0"/> ▲ ▼		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended				
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		16 ▲ ▼ ■ 33 ▲ ▼			
		Photometric <input type="text" value="1 Wavelength"/> ▼			
		Wavelength		Primary	<input type="text" value="570"/> ▼
				Secondary	<input type="text" value=""/> ▼

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Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: 1- Point ▼ Repls: 3 ▲▼ Curve Type: Linear ▼ Stability: 0 ▲▼ days				
Calibrator Selection N° 1 ▲▼ 1 # ▲▼ Calibrator: CAL Ferro Conc: @			R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Repls Range: 100 % Min Cal Resp: 0 mAbs ■ Reagent Blank Factor Change: 100 % ■ Auto Reagent Blank by Bottle M-Point Curve Fit: 0 %	

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="FOSFATASE"/>	Test Code: <input type="text" value="FAL79"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="U/L"/>	Decimal Points			
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="text" value="0"/>			
	<input type="radio"/> Extended				
		Methodology			
		Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="21"/>	<input type="text" value="25"/>		
		Photometric <input type="text" value="2 Wavelength"/>			
		Wavelength	Primary	<input type="text" value="405"/>	<input type="button" value="▼"/>
			Secondary	<input type="text" value="700"/>	<input type="button" value="▼"/>

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 FAL1	160	0	15 tests	0 days
R2 FAL2	40	0	15 tests	0 days

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FOSFATASE ALCALINA Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="27"/></td> <td><input type="text" value="100"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="27"/></td> <td><input type="text" value="100"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="27"/></td> <td><input type="text" value="100"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="27"/>	<input type="text" value="100"/>	Female:	<input type="text" value="27"/>	<input type="text" value="100"/>	Other:	<input type="text" value="27"/>	<input type="text" value="100"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="27"/>	<input type="text" value="100"/>														
Female:	<input type="text" value="27"/>	<input type="text" value="100"/>														
Other:	<input type="text" value="27"/>	<input type="text" value="100"/>														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Icterus</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Lipemia</td> <td><input type="text" value="0"/></td> </tr> </table>		Hemolysis	<input type="text" value="0"/>	Icterus	<input type="text" value="0"/>	Lipemia	<input type="text" value="0"/>	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>						
Hemolysis	<input type="text" value="0"/>															
Icterus	<input type="text" value="0"/>															
Lipemia	<input type="text" value="0"/>															
Valid Range: <input type="text" value="0"/> <input type="text" value="1500"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="40"/> %												
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off																

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: <input type="text" value="M-Point"/> ▼ Repts: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days													
Calibrator Selection <table border="1"> <tr> <td>N°</td> <td>Calibrator</td> <td>Conc</td> </tr> <tr> <td>2 ▲▼</td> <td>1 # ▲▼</td> <td>Calibra 1H @₁</td> </tr> <tr> <td>2</td> <td># ▲▼</td> <td>Calibra 2H @₂</td> </tr> </table>			N°	Calibrator	Conc	2 ▲▼	1 # ▲▼	Calibra 1H @ ₁	2	# ▲▼	Calibra 2H @ ₂	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Reps Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc											
2 ▲▼	1 # ▲▼	Calibra 1H @ ₁											
2	# ▲▼	Calibra 2H @ ₂											
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle										

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="FÓSFORO UV"/>	Test Code: <input type="text" value="FOS12"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/> ▼				
Reporting Unit:	<input type="text" value="mg/dL"/>	Decimal Points	<input type="text" value="2"/> ▲ ▼		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended				
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="24"/> ▲ ▼	<input type="checkbox"/>	<input type="text" value="0"/> ▲ ▼	
		Photometric <input type="text" value="2 Wavelength"/> ▼			
		Wavelength	Primary	<input type="text" value="340"/> ▼	
			Secondary	<input type="text" value="660"/> ▼	

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FÓSFORO UV Liquiform

Measure	Sample	Reagent	Ranges	Calibration																																																																											
<div style="display: flex; justify-content: space-between;"> <div style="width: 48%;"> <p>Sample Volumes</p> <table style="width:100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">Sample Vol.</th> <th style="text-align: center;">Dilution Sample Vol.</th> <th style="text-align: center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">2 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> </tbody> </table> <p>Diluent Code <input style="width: 100px;" type="text" value="Water"/></p> <p>Diluent Warming Limit <input style="width: 50px;" type="text" value="0"/> Tests <input type="text"/></p> </div> <div style="width: 48%;"> <p>Condition</p> <table style="width:100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">1</th> <th style="text-align: center;">2</th> <th style="text-align: center;">3</th> <th style="text-align: center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> </tbody> </table> </div> </div>						Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	2 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	3	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	4	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>		1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	2 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
3	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
4	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 FOS	200 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/> tests	0 <input type="text"/> days
		<input checked="" type="checkbox"/> Stirring		
		<input type="checkbox"/> Stirring		

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FÓSFORO UV Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="2.5"/></td> <td><input type="text" value="4.8"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="2.5"/></td> <td><input type="text" value="4.8"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="2.5"/></td> <td><input type="text" value="4.8"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="2.5"/>	<input type="text" value="4.8"/>	Female:	<input type="text" value="2.5"/>	<input type="text" value="4.8"/>	Other:	<input type="text" value="2.5"/>	<input type="text" value="4.8"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="2.5"/>	<input type="text" value="4.8"/>														
Female:	<input type="text" value="2.5"/>	<input type="text" value="4.8"/>														
Other:	<input type="text" value="2.5"/>	<input type="text" value="4.8"/>														
Valid Range: <input type="text" value="0"/> <input type="text" value="20"/>		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>												
Serum Index Limit Hemolysis: <input type="text" value="0"/> Icterus: <input type="text" value="0"/> Lipemia: <input type="text" value="0"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: <input type="text" value="M-Point"/> ▼ Repls: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days													
Calibrator Selection <table border="1"> <thead> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td><input type="text" value="2"/> ▲▼</td> <td><input type="text" value="1"/> # <input type="text" value="Calibra 2H"/></td> <td><input type="text" value="@1"/></td> </tr> <tr> <td><input type="text" value="2"/> ▲▼</td> <td><input type="text" value="2"/> # <input type="text" value="Calibra 1H"/></td> <td><input type="text" value="@2"/></td> </tr> </tbody> </table>			N°	Calibrator	Conc	<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # <input type="text" value="Calibra 2H"/>	<input type="text" value="@1"/>	<input type="text" value="2"/> ▲▼	<input type="text" value="2"/> # <input type="text" value="Calibra 1H"/>	<input type="text" value="@2"/>	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repls Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc											
<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # <input type="text" value="Calibra 2H"/>	<input type="text" value="@1"/>											
<input type="text" value="2"/> ▲▼	<input type="text" value="2"/> # <input type="text" value="Calibra 1H"/>	<input type="text" value="@2"/>											
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle										

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="FRUTOSAMINA"/>	Test Code: <input type="text" value="FRU97"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>			
Reporting Unit:	<input type="text" value="µmol/L"/>			
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended			
Decimal Points	<input type="text" value="0"/>			
Methodology				
Type:		<input checked="" type="radio"/> End Point <input type="radio"/> Rate		
Measuring Point:				
23		33		
Photometric		<input type="text" value="2 Wavelength"/>		
Wavelength		Primary <input type="text" value="546"/>		
		Secondary <input type="text" value="700"/>		

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Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: 1- Point Repls: 3 Curve Type: Linear Stability: 0 days				
Calibrator Selection N°: 1 Calibrator: CAL Fru Conc: @			R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Repls Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit: 0 %	
<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle				

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="Gama GT"/>	Test Code: <input type="text" value="GGT105"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology		
Reporting Unit:	<input type="text" value="U/L"/>	Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:		
Decimal Points	<input type="text" value="2"/>	<input type="text" value="21"/>	<input type="text" value="25"/>	
	<input type="button" value="▲"/> <input type="button" value="▼"/>	<input type="button" value="▲"/> <input type="button" value="▼"/>	<input type="button" value="▲"/> <input type="button" value="▼"/>	
		Photometric	<input type="text" value="2 Wavelength"/>	
		Wavelength	Primary	<input type="text" value="405"/>
			Secondary	<input type="text" value="700"/>

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GAMA GT Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="7"/></td> <td><input type="text" value="58"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="5"/></td> <td><input type="text" value="39"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="7"/></td> <td><input type="text" value="58"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="7"/>	<input type="text" value="58"/>	Female:	<input type="text" value="5"/>	<input type="text" value="39"/>	Other:	<input type="text" value="7"/>	<input type="text" value="58"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="7"/>	<input type="text" value="58"/>														
Female:	<input type="text" value="5"/>	<input type="text" value="39"/>														
Other:	<input type="text" value="7"/>	<input type="text" value="58"/>														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Icterus</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Lipemia</td> <td><input type="text" value="0"/></td> </tr> </table>		Hemolysis	<input type="text" value="0"/>	Icterus	<input type="text" value="0"/>	Lipemia	<input type="text" value="0"/>	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>						
Hemolysis	<input type="text" value="0"/>															
Icterus	<input type="text" value="0"/>															
Lipemia	<input type="text" value="0"/>															
Valid Range: <input type="text" value="0"/> <input type="text" value="700"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="40"/> %												
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off																

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: <input type="text" value="M-Point"/> ▼ Repeats: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days													
Calibrator Selection <table border="1"> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> <tr> <td><input type="text" value="2"/> ▲▼</td> <td><input type="text" value="1"/> # ▲▼</td> <td><input type="text" value="Calibra 1H"/> @₁</td> </tr> <tr> <td></td> <td><input type="text" value="2"/> # ▲▼</td> <td><input type="text" value="Calibra 2H"/> @₂</td> </tr> </table>		N°	Calibrator	Conc	<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # ▲▼	<input type="text" value="Calibra 1H"/> @ ₁		<input type="text" value="2"/> # ▲▼	<input type="text" value="Calibra 2H"/> @ ₂	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repts Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %		
N°	Calibrator	Conc											
<input type="text" value="2"/> ▲▼	<input type="text" value="1"/> # ▲▼	<input type="text" value="Calibra 1H"/> @ ₁											
	<input type="text" value="2"/> # ▲▼	<input type="text" value="Calibra 2H"/> @ ₂											
<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle													

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Test N°	<input type="text" value="#"/>				<input type="button" value="OK"/>
Test Name:	<input type="text" value="GLICOSE HK"/>	Test Code:	<input type="text" value="GHK85"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration	
Sample Type:	<input type="text" value="Serum"/>	Methodology			
Reporting Unit:	<input type="text" value="mg/dL"/>	Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:			
Decimal Points	<input type="text" value="0"/>	16 <input type="text"/> 23			
		Photometric <input type="text" value="2 Wavelength"/>			
		Wavelength Primary <input type="text" value="340"/>			
		Secondary <input type="text" value="375"/>			

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 GHK1	240	0	15 tests	0 days
R2 GHK2	60	0	15 tests	0 days
		■ Stirring		

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="GLICOSE PAP"/>	Test Code: <input type="text" value="GLI84"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology		
Reporting Unit:	<input type="text" value="mg/dL"/>	Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:		
Decimal Points	<input type="text" value="0"/>	2 <input type="text" value="33"/>		
		Photometric <input type="text" value="2 Wavelength"/>		
		Wavelength		
		Primary <input type="text" value="510"/>		
		Secondary <input type="text" value="660"/>		

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 GLI	200 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/> tests	0 <input type="text"/> days
		<input checked="" type="checkbox"/> Stirring		
		<input type="checkbox"/> Stirring		

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GLICOSE PAP Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table style="width:100%; border-collapse: collapse;"> <tr> <td></td> <td style="text-align:center;">Lower</td> <td style="text-align:center;">Upper</td> </tr> <tr> <td>Male:</td> <td style="text-align:center;"><input type="text" value="65"/></td> <td style="text-align:center;"><input type="text" value="99"/></td> </tr> <tr> <td>Female:</td> <td style="text-align:center;"><input type="text" value="65"/></td> <td style="text-align:center;"><input type="text" value="99"/></td> </tr> <tr> <td>Other:</td> <td style="text-align:center;"><input type="text" value="65"/></td> <td style="text-align:center;"><input type="text" value="99"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="65"/>	<input type="text" value="99"/>	Female:	<input type="text" value="65"/>	<input type="text" value="99"/>	Other:	<input type="text" value="65"/>	<input type="text" value="99"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="65"/>	<input type="text" value="99"/>														
Female:	<input type="text" value="65"/>	<input type="text" value="99"/>														
Other:	<input type="text" value="65"/>	<input type="text" value="99"/>														
		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs	Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>													
Valid Range: <input type="text" value="0"/> <input type="text" value="500"/>		Non-Linear Limit: <input type="text" value="0"/> %														
Serum Index Limit Hemolysis: <input type="text" value="0"/> Icterus: <input type="text" value="0"/> Lipemia: <input type="text" value="0"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration																		
Calibration Method: <input type="text" value="M-Point"/> ▼ Repeats: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days																						
Calibrator Selection <table style="width:100%; border-collapse: collapse;"> <tr> <td>N°</td> <td>Calibrator</td> <td>Conc</td> </tr> <tr> <td style="text-align:center;">2</td> <td style="text-align:center;">#</td> <td style="text-align:center;">#</td> </tr> <tr> <td style="text-align:center;">1</td> <td style="text-align:center;">#</td> <td style="text-align:center;">#</td> </tr> <tr> <td style="text-align:center;">2</td> <td style="text-align:center;">#</td> <td style="text-align:center;">#</td> </tr> <tr> <td></td> <td>Calibra 2H</td> <td>@1</td> </tr> <tr> <td></td> <td>Calibra 1H</td> <td>@2</td> </tr> </table>			N°	Calibrator	Conc	2	#	#	1	#	#	2	#	#		Calibra 2H	@1		Calibra 1H	@2	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repts Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc																				
2	#	#																				
1	#	#																				
2	#	#																				
	Calibra 2H	@1																				
	Calibra 1H	@2																				
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle																			

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="HDL LE"/>	Test Code: <input type="text" value="HDL98"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="mg/dL"/>	Decimal Points <input type="text" value="0"/>	Methodology		
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="radio"/> Extended	Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate		
			Measuring Point:		
			<input type="text" value="16"/>	<input type="text" value="33"/>	
			Photometric <input type="text" value="2 Wavelength"/>		
			Wavelength	Primary	<input type="text" value="600"/>
				Secondary	<input type="text" value="700"/>

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HDL LE

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="35"/></td> <td><input type="text" value="65"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="35"/></td> <td><input type="text" value="65"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="35"/></td> <td><input type="text" value="65"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="35"/>	<input type="text" value="65"/>	Female:	<input type="text" value="35"/>	<input type="text" value="65"/>	Other:	<input type="text" value="35"/>	<input type="text" value="65"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="35"/>	<input type="text" value="65"/>														
Female:	<input type="text" value="35"/>	<input type="text" value="65"/>														
Other:	<input type="text" value="35"/>	<input type="text" value="65"/>														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Icterus</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Lipemia</td> <td><input type="text" value="0"/></td> </tr> </table>		Hemolysis	<input type="text" value="0"/>	Icterus	<input type="text" value="0"/>	Lipemia	<input type="text" value="0"/>	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>						
Hemolysis	<input type="text" value="0"/>															
Icterus	<input type="text" value="0"/>															
Lipemia	<input type="text" value="0"/>															
Valid Range: <input type="text" value="2.5"/> <input type="text" value="200"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> %												
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off																

Measure	Sample	Reagent	Ranges	Calibration						
Calibration Method: <input type="text" value="1- Point"/> ▼ Repeats: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days										
Calibrator Selection <table border="1"> <tr> <td>N°</td> <td>Calibrator</td> <td>Conc</td> </tr> <tr> <td><input type="text" value="1"/> ▲▼</td> <td><input type="text" value="1"/> #</td> <td><input type="text" value="CAL HDL"/> @</td> </tr> </table>			N°	Calibrator	Conc	<input type="text" value="1"/> ▲▼	<input type="text" value="1"/> #	<input type="text" value="CAL HDL"/> @	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repeats Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc								
<input type="text" value="1"/> ▲▼	<input type="text" value="1"/> #	<input type="text" value="CAL HDL"/> @								
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle							

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="CAP LIG FER"/>	Test Code: <input type="text" value="IBC92"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="µg/dL"/>	Decimal Points	<input type="text" value="0"/>		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended				
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="16"/>	■	<input type="text" value="33"/>	
		Photometric <input type="text" value="1 Wavelength"/>			
		Wavelength	Primary	<input type="text" value="570"/>	
			Secondary	<input type="text"/>	

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 IBC1 <input type="text"/>	190 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
<input type="checkbox"/> Stirring				
R2 IBC2 <input type="text"/>	50 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
<input type="checkbox"/> Stirring				

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IBC Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="250"/></td> <td><input type="text" value="450"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="250"/></td> <td><input type="text" value="450"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="250"/></td> <td><input type="text" value="450"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="250"/>	<input type="text" value="450"/>	Female:	<input type="text" value="250"/>	<input type="text" value="450"/>	Other:	<input type="text" value="250"/>	<input type="text" value="450"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="250"/>	<input type="text" value="450"/>														
Female:	<input type="text" value="250"/>	<input type="text" value="450"/>														
Other:	<input type="text" value="250"/>	<input type="text" value="450"/>														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td><input type="text"/></td> </tr> <tr> <td>Icterus</td> <td><input type="text"/></td> </tr> <tr> <td>Lipemia</td> <td><input type="text"/></td> </tr> </table>		Hemolysis	<input type="text"/>	Icterus	<input type="text"/>	Lipemia	<input type="text"/>	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>						
Hemolysis	<input type="text"/>															
Icterus	<input type="text"/>															
Lipemia	<input type="text"/>															
Valid Range: <input type="text" value="0"/> <input type="text" value="500"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> %												
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off																

Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: <input type="text" value="1- Point"/> ▼ Repls: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days				
Calibrator Selection N° <input type="text" value="1"/> ▲▼ <input type="text" value="1"/> # <input type="text" value="1"/> ▲▼ Calibrator: <input type="text" value="CAL Ferro"/> Conc: <input type="text" value="@"/>			R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repls Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle	

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="LDH"/>	Test Code: <input type="text" value="LDH86"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology			
Reporting Unit:	<input type="text" value="U/L"/>	Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate			
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:			
	Decimal Points: <input type="text" value="0"/>	21 <input type="text" value="▲▼"/> 25 <input type="text" value="▲▼"/>			
		Photometric: <input type="text" value="2 Wavelength"/>			
		Wavelength: Primary <input type="text" value="340"/> <input type="text" value="▼"/>			
		Secondary <input type="text" value="700"/> <input type="text" value="▼"/>			

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Diluent Vol.	Warming Limit	Stability
R1 LDH1 <input type="button" value="▲"/> <input type="button" value="▼"/>	160 <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="button" value="▲"/> <input type="button" value="▼"/>	<input type="checkbox"/> Stirring	15 <input type="button" value="▲"/> <input type="button" value="▼"/> tests
R2 LDH2 <input type="button" value="▲"/> <input type="button" value="▼"/>	40 <input type="button" value="▲"/> <input type="button" value="▼"/>	0 <input type="button" value="▲"/> <input type="button" value="▼"/>	<input type="checkbox"/> Stirring	15 <input type="button" value="▲"/> <input type="button" value="▼"/> tests
				0 <input type="button" value="▲"/> <input type="button" value="▼"/> days
				0 <input type="button" value="▲"/> <input type="button" value="▼"/> days

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LDH Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="200"/></td> <td><input type="text" value="480"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="200"/></td> <td><input type="text" value="480"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="200"/></td> <td><input type="text" value="480"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="200"/>	<input type="text" value="480"/>	Female:	<input type="text" value="200"/>	<input type="text" value="480"/>	Other:	<input type="text" value="200"/>	<input type="text" value="480"/>	Reaction Slope <input checked="" type="radio"/> Negative <input type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="200"/>	<input type="text" value="480"/>														
Female:	<input type="text" value="200"/>	<input type="text" value="480"/>														
Other:	<input type="text" value="200"/>	<input type="text" value="480"/>														
Valid Range: <input type="text" value="0"/> <input type="text" value="2000"/>		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>												
Serum Index Limit Hemolysis: <input type="text" value="0"/> Icterus: <input type="text" value="0"/> Lipemia: <input type="text" value="0"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="40"/> % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration															
Calibration Method: <input type="text" value="M-Point"/> ▼ Repeats: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days																			
Calibrator Selection <table border="1"> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> <tr> <td>2</td> <td># <input type="text" value="1"/></td> <td><input type="text" value="@1"/></td> </tr> <tr> <td></td> <td>Calibra 1H</td> <td></td> </tr> <tr> <td>2</td> <td># <input type="text" value="2"/></td> <td><input type="text" value="@2"/></td> </tr> <tr> <td></td> <td>Calibra 2H</td> <td></td> </tr> </table>			N°	Calibrator	Conc	2	# <input type="text" value="1"/>	<input type="text" value="@1"/>		Calibra 1H		2	# <input type="text" value="2"/>	<input type="text" value="@2"/>		Calibra 2H		R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="-100"/> mAbs Cal Repts Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc																	
2	# <input type="text" value="1"/>	<input type="text" value="@1"/>																	
	Calibra 1H																		
2	# <input type="text" value="2"/>	<input type="text" value="@2"/>																	
	Calibra 2H																		
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle																

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="LDL"/>	Test Code: <input type="text" value="LDL111"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="mg/dL"/>	Decimal Points			
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="text" value="0"/>			
	<input type="radio"/> Extended	<input type="button" value="▲"/>			
		<input type="button" value="▼"/>			
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="16"/>	<input type="text" value="33"/>		
		<input type="button" value="▲"/>	<input type="button" value="▲"/>		
		<input type="button" value="▼"/>	<input type="button" value="▼"/>		
		Photometric	<input type="text" value="2 Wavelength"/>		
		Wavelength	Primary	<input type="text" value="546"/>	<input type="button" value="▼"/>
			Secondary	<input type="text" value="700"/>	<input type="button" value="▼"/>

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LDL Liquiform

Measure	Sample	Reagent	Ranges	Calibration																																																																											
<div style="display: flex; justify-content: space-between;"> <div style="width: 48%;"> <p>Sample Volumes</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">Sample Vol.</th> <th style="text-align: center;">Dilution Sample Vol.</th> <th style="text-align: center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">3 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> </tbody> </table> <p>Diluent Code <input style="width: 100px;" type="text" value="Water"/></p> <p>Diluent Warming Limit <input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/> Tests</p> </div> <div style="width: 48%;"> <p>Condition</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">1</th> <th style="text-align: center;">2</th> <th style="text-align: center;">3</th> <th style="text-align: center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> </tbody> </table> </div> </div>						Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	3 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	2	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	3	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	4	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>		1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	3 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
2	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
3	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
4	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration
Reagent Code Volume Diluente Vol. Warming Limit Stability				
R1 LDL1 <input type="text" value="▲"/> <input type="text" value="▼"/> 225 <input type="text" value="▲"/> <input type="text" value="▼"/> 0 <input type="text" value="▲"/> <input type="text" value="▼"/> <input checked="" type="checkbox"/> Stirring 15 <input type="text" value="▲"/> <input type="text" value="▼"/> tests 0 <input type="text" value="▲"/> <input type="text" value="▼"/> days				
R2 LDL2 <input type="text" value="▲"/> <input type="text" value="▼"/> 75 <input type="text" value="▲"/> <input type="text" value="▼"/> 0 <input type="text" value="▲"/> <input type="text" value="▼"/> <input checked="" type="checkbox"/> Stirring 15 <input type="text" value="▲"/> <input type="text" value="▼"/> tests 0 <input type="text" value="▲"/> <input type="text" value="▼"/> days				

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LDL Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="7"/></td> <td><input type="text" value="129"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="7"/></td> <td><input type="text" value="129"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="7"/></td> <td><input type="text" value="129"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="7"/>	<input type="text" value="129"/>	Female:	<input type="text" value="7"/>	<input type="text" value="129"/>	Other:	<input type="text" value="7"/>	<input type="text" value="129"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="7"/>	<input type="text" value="129"/>														
Female:	<input type="text" value="7"/>	<input type="text" value="129"/>														
Other:	<input type="text" value="7"/>	<input type="text" value="129"/>														
Valid Range: <input type="text" value="6.6"/> <input type="text" value="992"/>		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept <input type="text" value="0"/>												
Serum Index Limit Hemolysis <input type="text" value="0"/> Icterus <input type="text" value="0"/> Lipemia <input type="text" value="0"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit <input type="text" value="0"/> Equation <input type="text" value="none"/> ▼ Judge Point <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: <input type="text" value="1- Point"/> ▼ Repls: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability <input type="text" value="0"/> ▲▼ days				
Calibrator Selection N° <input type="text" value="1"/> ▲▼ <input type="text" value="1"/> # <input type="text" value="CAL LDL"/> <input type="text" value="@"/>			R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repls Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit <input type="text" value="0"/> %	
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle	

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="LIPASE"/>	Test Code: <input type="text" value="LIP107"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/> ▼				
Reporting Unit:	<input type="text" value="U/L"/>	Decimal Points			
		<input type="text" value="0"/> ▲ ▼			
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended				
		Methodology			
		Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="23"/> ▲ ▼	<input type="text" value="33"/> ▲ ▼		
		Photometric <input type="text" value="2 Wavelength"/> ▼			
		Wavelength		Primary	<input type="text" value="570"/> ▼
				Secondary	<input type="text" value="700"/> ▼

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Lipase Liquiform

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="13"/></td> <td><input type="text" value="60"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="13"/></td> <td><input type="text" value="60"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="13"/></td> <td><input type="text" value="60"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="13"/>	<input type="text" value="60"/>	Female:	<input type="text" value="13"/>	<input type="text" value="60"/>	Other:	<input type="text" value="13"/>	<input type="text" value="60"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	<input type="text" value="13"/>	<input type="text" value="60"/>														
Female:	<input type="text" value="13"/>	<input type="text" value="60"/>														
Other:	<input type="text" value="13"/>	<input type="text" value="60"/>														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Icterus</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Lipemia</td> <td><input type="text" value="0"/></td> </tr> </table>		Hemolysis	<input type="text" value="0"/>	Icterus	<input type="text" value="0"/>	Lipemia	<input type="text" value="0"/>	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>						
Hemolysis	<input type="text" value="0"/>															
Icterus	<input type="text" value="0"/>															
Lipemia	<input type="text" value="0"/>															
Valid Range: <input type="text" value="3"/> <input type="text" value="300"/>		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="40"/> %												
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off																

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: <input type="text" value="M-Point"/> ▼ Repls: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days													
Calibrator Selection <table border="1"> <thead> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td><input type="text" value="2"/> ▲▼</td> <td><input type="text" value="1 #"/> ▲▼</td> <td><input type="text" value="Calibra 1H"/> @₁</td> </tr> <tr> <td></td> <td><input type="text" value="2 #"/> ▲▼</td> <td><input type="text" value="Calibra 2H"/> @₂</td> </tr> </tbody> </table>			N°	Calibrator	Conc	<input type="text" value="2"/> ▲▼	<input type="text" value="1 #"/> ▲▼	<input type="text" value="Calibra 1H"/> @ ₁		<input type="text" value="2 #"/> ▲▼	<input type="text" value="Calibra 2H"/> @ ₂	R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repls Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %	
N°	Calibrator	Conc											
<input type="text" value="2"/> ▲▼	<input type="text" value="1 #"/> ▲▼	<input type="text" value="Calibra 1H"/> @ ₁											
	<input type="text" value="2 #"/> ▲▼	<input type="text" value="Calibra 2H"/> @ ₂											
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle										

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="MAGNÉSIO"/>	Test Code: <input type="text" value="Mg50"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>	Methodology		
Reporting Unit:	<input type="text" value="mg/dL"/>	Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate		
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended	Measuring Point:		
Decimal Points	<input type="text" value="2"/>	<input type="text" value="11"/>		
	<input type="button" value="▲"/> <input type="button" value="▼"/>	<input type="text" value="0"/>		
		<input type="button" value="▲"/> <input type="button" value="▼"/>		
		Photometric <input type="text" value="1 Wavelength"/>		
		Wavelength Primary <input type="text" value="510"/>		
		Secondary <input type="text" value=""/>		

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MAGNÉSIO

Measure	Sample	Reagent	Ranges	Calibration												
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td>1.9</td> <td>2.5</td> </tr> <tr> <td>Female:</td> <td>1.9</td> <td>2.5</td> </tr> <tr> <td>Other:</td> <td>1.9</td> <td>2.5</td> </tr> </table>			Lower	Upper	Male:	1.9	2.5	Female:	1.9	2.5	Other:	1.9	2.5	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper														
Male:	1.9	2.5														
Female:	1.9	2.5														
Other:	1.9	2.5														
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td>0</td> </tr> <tr> <td>Icterus</td> <td>0</td> </tr> <tr> <td>Lipemia</td> <td>0</td> </tr> </table>		Hemolysis	0	Icterus	0	Lipemia	0	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 3500 mAbs		Correction Constant Slope: 1 Intercept: 0						
Hemolysis	0															
Icterus	0															
Lipemia	0															
Valid Range: 0 4.5		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: 0 Equation: none Judge Point: 0		Non-Linear Limit: 0 % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off												

Measure	Sample	Reagent	Ranges	Calibration									
Calibration Method: M-Point Repeats: 3 Curve Type: Linear Stability: 0 days													
Calibrator Selection <table border="1"> <tr> <th>N°</th> <th>Calibrator</th> <th>Conc</th> </tr> <tr> <td>2</td> <td>Calibra 2H</td> <td>@1</td> </tr> <tr> <td>2</td> <td>Calibra 1H</td> <td>@2</td> </tr> </table>			N°	Calibrator	Conc	2	Calibra 2H	@1	2	Calibra 1H	@2	R-Blank + Calibration Limits R-Blank Limit: 3500 mAbs Cal Repts Range: 100 % Min Cal Resp: 0 mAbs Factor Change: 100 % M-Point Curve Fit: 0 %	
N°	Calibrator	Conc											
2	Calibra 2H	@1											
2	Calibra 1H	@2											
			<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle										

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="PROTEÍNAS"/>	Test Code: <input type="text" value="PRO99"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>			
Reporting Unit:	<input type="text" value="g/dL"/>	Decimal Points		
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="radio"/> Extended		
		Methodology		
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate		
		Measuring Point:		
		<input type="text" value="33"/>	<input type="text" value="0"/>	
		Photometric <input type="text" value="2 Wavelength"/>		
		Wavelength	Primary	<input type="text" value="546"/>
			Secondary	<input type="text" value="700"/>

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PROTEÍNAS TOTAIS

Measure	Sample	Reagent	Ranges	Calibration
Normal Range Lower Upper Male: <input type="text" value="6"/> <input type="text" value="8"/> Female: <input type="text" value="6"/> <input type="text" value="8"/> Other: <input type="text" value="6"/> <input type="text" value="8"/> Valid Range: <input type="text" value="1"/> <input type="text" value="14"/>		Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
Serum Index Limit Hemolysis <input type="text" value="0"/> Icterus <input type="text" value="0"/> Lipemia <input type="text" value="0"/>		Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept <input type="text" value="0"/>
		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit <input type="text" value="0"/> Equation <input type="text" value="none"/> ▼ Judge Point <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off

Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: <input type="text" value="M-Point"/> ▼ Repts: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability <input type="text" value="0"/> ▲▼ days				
Calibrator Selection N° Calibrator Conc 2 <input type="text" value="1"/> # <input type="text" value="Calibra 2H"/> <input type="text" value="@1"/> 2 <input type="text" value="2"/> # <input type="text" value="Calibra 1H"/> <input type="text" value="@2"/>		R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Reps Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit <input type="text" value="0"/> %		
		<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle		

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="SENSIPROT"/>	Test Code: <input type="text" value="SEN36"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="mg/dL"/>	Decimal Points			
		<input type="text" value="0"/>			
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="radio"/> Extended			
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="33"/>	<input type="checkbox"/>	<input type="text" value="0"/>	
		Photometric <input type="text" value="2 Wavelength"/>			
		Wavelength	Primary	<input type="text" value="600"/>	<input type="button" value="▼"/>
			Secondary	<input type="text" value="660"/>	<input type="button" value="▼"/>

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SENSIPROT

Measure	Sample	Reagent	Ranges	Calibration																																																																											
<p>Sample Volumes</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">Sample Vol.</th> <th style="text-align: center;">Dilution Sample Vol.</th> <th style="text-align: center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">9 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> <td style="text-align: center;">0 <input type="text" value="▲"/> <input type="text" value="▼"/></td> </tr> </tbody> </table> <p>Diluent Code <input style="width: 100px;" type="text" value="Water"/></p> <p>Diluent Warming Limit <input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/> Tests</p>			Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	9 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	2	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	3	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	4	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	<p>Condition</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">1</th> <th style="text-align: center;">2</th> <th style="text-align: center;">3</th> <th style="text-align: center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> </tbody> </table>				1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	9 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
2	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
3	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
4	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>	0 <input type="text" value="▲"/> <input type="text" value="▼"/>																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration	
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability	
R1 SEN	<input style="width: 50px;" type="text" value="180"/> <input type="text" value="▲"/> <input type="text" value="▼"/>	<input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/>	<input checked="" type="checkbox"/> Stirring	<input style="width: 50px;" type="text" value="15"/> <input type="text" value="▲"/> <input type="text" value="▼"/> tests	<input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/> days
R2	<input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/>	<input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/>	<input type="checkbox"/> Stirring	<input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/> tests	<input style="width: 50px;" type="text" value="0"/> <input type="text" value="▲"/> <input type="text" value="▼"/> days

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SENSIPROT

Measure	Sample	Reagent	Ranges	Calibration															
Normal Range <table border="1"> <tr> <td></td> <td>Lower</td> <td>Upper</td> </tr> <tr> <td>Male:</td> <td><input type="text" value="1"/></td> <td><input type="text" value="15"/></td> </tr> <tr> <td>Female:</td> <td><input type="text" value="1"/></td> <td><input type="text" value="15"/></td> </tr> <tr> <td>Other:</td> <td><input type="text" value="1"/></td> <td><input type="text" value="15"/></td> </tr> <tr> <td>Valid Range:</td> <td><input type="text" value="0"/></td> <td><input type="text" value="100"/></td> </tr> </table>			Lower	Upper	Male:	<input type="text" value="1"/>	<input type="text" value="15"/>	Female:	<input type="text" value="1"/>	<input type="text" value="15"/>	Other:	<input type="text" value="1"/>	<input type="text" value="15"/>	Valid Range:	<input type="text" value="0"/>	<input type="text" value="100"/>	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive		
	Lower	Upper																	
Male:	<input type="text" value="1"/>	<input type="text" value="15"/>																	
Female:	<input type="text" value="1"/>	<input type="text" value="15"/>																	
Other:	<input type="text" value="1"/>	<input type="text" value="15"/>																	
Valid Range:	<input type="text" value="0"/>	<input type="text" value="100"/>																	
Serum Index Limit <table border="1"> <tr> <td>Hemolysis</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Icterus</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Lipemia</td> <td><input type="text" value="0"/></td> </tr> </table>		Hemolysis	<input type="text" value="0"/>	Icterus	<input type="text" value="0"/>	Lipemia	<input type="text" value="0"/>	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="3500"/> mAbs		Correction Constant Slope: <input type="text" value="1"/> Intercept: <input type="text" value="0"/>									
Hemolysis	<input type="text" value="0"/>																		
Icterus	<input type="text" value="0"/>																		
Lipemia	<input type="text" value="0"/>																		
		Prozone Limit <input checked="" type="radio"/> Above <input type="radio"/> Below Limit: <input type="text" value="0"/> Equation: <input type="text" value="none"/> ▼ Judge Point: <input type="text" value="0"/> ▲▼		Non-Linear Limit: <input type="text" value="0"/> % Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off															

Measure	Sample	Reagent	Ranges	Calibration
Calibration Method: <input type="text" value="1- Point"/> ▼ Repls: <input type="text" value="3"/> ▲▼ Curve Type: <input type="text" value="Linear"/> ▼ Stability: <input type="text" value="0"/> ▲▼ days				
Calibrator Selection N° <input type="text" value="1"/> ▲▼ <input type="text" value="1"/> # <input type="text" value=""/> ▲▼ <input type="text" value="CAL Sens"/> <input type="text" value="@"/>		R-Blank + Calibration Limits R-Blank Limit: <input type="text" value="3500"/> mAbs Cal Repls Range: <input type="text" value="100"/> % Min Cal Resp: <input type="text" value="0"/> mAbs Factor Change: <input type="text" value="100"/> % M-Point Curve Fit: <input type="text" value="0"/> %		
		<input checked="" type="checkbox"/> Reagent Blank <input checked="" type="checkbox"/> Auto Reagent Blank by Bottle		

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>			
Test Name:	<input type="text" value="TRIGLICÉRIDE"/>	Test Code: <input type="text" value="TRI87"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>	
Measure		Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>				
Reporting Unit:	<input type="text" value="mg/dL"/>	Decimal Points			
Reaction Cycle	<input checked="" type="radio"/> Standard	<input type="text" value="0"/>			
	<input type="radio"/> Extended				
		Methodology			
		Type: <input checked="" type="radio"/> End Point <input type="radio"/> Rate			
		Measuring Point:			
		<input type="text" value="33"/>	<input type="checkbox"/>	<input type="text" value="0"/>	
		Photometric <input type="text" value="2 Wavelength"/>			
		Wavelength	Primary	<input type="text" value="510"/>	<input type="button" value="▼"/>
			Secondary	<input type="text" value="660"/>	<input type="button" value="▼"/>

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TRIGLICÉRIDES Liquiform

Measure	Sample	Reagent	Ranges	Calibration																																																																											
<div style="display: flex; justify-content: space-between;"> <div style="width: 48%;"> <p>Sample Volumes</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">Sample Vol.</th> <th style="text-align: center;">Dilution Sample Vol.</th> <th style="text-align: center;">Dilution Diluent Vol.</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">2 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> </tr> </tbody> </table> <p style="margin-top: 10px;">Diluent Code Water</p> <p>Diluent Warming Limit 0 <input type="text"/> Tests <input type="text"/></p> </div> <div style="width: 48%;"> <p>Condition</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">1</th> <th style="text-align: center;">2</th> <th style="text-align: center;">3</th> <th style="text-align: center;">4</th> </tr> </thead> <tbody> <tr> <td>First Run</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Below N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Above N-Range</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> Panic L</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Panic H</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Noise</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> Prozone</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input type="checkbox"/> HIGH!</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> <tr> <td><input type="checkbox"/> ABS!</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Samp. Vol Reduction</td> <td style="text-align: center;">●</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> <td style="text-align: center;">○</td> </tr> </tbody> </table> </div> </div>						Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.	1	2 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	3	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>	4	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>		1	2	3	4	First Run	●	○	○	○	<input type="checkbox"/> Below N-Range					<input type="checkbox"/> Above N-Range					<input type="checkbox"/> Panic L	○	○	○	○	<input type="checkbox"/> Panic H	○	○	○	○	<input type="checkbox"/> Noise	○	○	○	○	<input type="checkbox"/> Prozone					<input type="checkbox"/> HIGH!	○	○	○	○	<input type="checkbox"/> ABS!					Samp. Vol Reduction	●	○	○	○
	Sample Vol.	Dilution Sample Vol.	Dilution Diluent Vol.																																																																												
1	2 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
3	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
4	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/>																																																																												
	1	2	3	4																																																																											
First Run	●	○	○	○																																																																											
<input type="checkbox"/> Below N-Range																																																																															
<input type="checkbox"/> Above N-Range																																																																															
<input type="checkbox"/> Panic L	○	○	○	○																																																																											
<input type="checkbox"/> Panic H	○	○	○	○																																																																											
<input type="checkbox"/> Noise	○	○	○	○																																																																											
<input type="checkbox"/> Prozone																																																																															
<input type="checkbox"/> HIGH!	○	○	○	○																																																																											
<input type="checkbox"/> ABS!																																																																															
Samp. Vol Reduction	●	○	○	○																																																																											

Measure	Sample	Reagent	Ranges	Calibration															
<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Reagent Code</th> <th style="width: 15%;">Volume</th> <th style="width: 15%;">Dilute Vol.</th> <th style="width: 15%;">Warming Limit</th> <th style="width: 15%;">Stability</th> </tr> </thead> <tbody> <tr> <td>R1 TRI</td> <td style="text-align: center;">200 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">15 <input type="text"/> tests</td> <td style="text-align: center;">0 <input type="text"/> days</td> </tr> <tr> <td>R2</td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/></td> <td style="text-align: center;">0 <input type="text"/> tests</td> <td style="text-align: center;">0 <input type="text"/> days</td> </tr> </tbody> </table> <p style="margin-top: 10px;"> <input checked="" type="checkbox"/> Stirring <input type="checkbox"/> Stirring </p>					Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability	R1 TRI	200 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days	R2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/> tests	0 <input type="text"/> days
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability															
R1 TRI	200 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days															
R2	0 <input type="text"/>	0 <input type="text"/>	0 <input type="text"/> tests	0 <input type="text"/> days															

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Test N°	<input type="text" value="#"/>	<input type="button" value="OK"/>		
Test Name:	<input type="text" value="URÉIA UV"/>	Test Code: <input type="text" value="UR104"/>	<input type="button" value="Print"/>	<input type="button" value="Cancel"/>
Measure	Sample	Reagent	Ranges	Calibration
Sample Type:	<input type="text" value="Serum"/>			
Reporting Unit:	<input type="text" value="mg/dL"/>			
Reaction Cycle	<input checked="" type="radio"/> Standard <input type="radio"/> Extended			
Decimal Points	<input type="text" value="0"/>			
		Methodology		
		Type: <input type="radio"/> End Point <input checked="" type="radio"/> Rate		
		Measuring Point:		
		<input type="text" value="20"/>	<input type="text" value="24"/>	
		Photometric <input type="text" value="2 Wavelength"/>		
		Wavelength	Primary	<input type="text" value="340"/>
			Secondary	<input type="text" value="700"/>

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Measure	Sample	Reagent	Ranges	Calibration
Reagent Code	Volume	Dilute Vol.	Warming Limit	Stability
R1 URÉIA1 <input type="text"/>	160 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
R2 URÉIA2 <input type="text"/>	40 <input type="text"/>	0 <input type="text"/>	15 <input type="text"/> tests	0 <input type="text"/> days
		■ Stirring		

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